

Product Datasheet

Recombinant Human RASSF1 His Protein NBP2-23397

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-23397

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-23397




NBP2-23397

Recombinant Human RASSF1 His Protein

| Product Information | |
|-----------------------------|--|
| Unit Size | 0.1 mg |
| Concentration | 1 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Preservative | No Preservative |
| Purity | >85%, by SDS-PAGE |
| Buffer | 20 mM Tris-HCl buffer (pH 8.0), 0.4M UREA, 10% glycerol |
| Target Molecular Weight | 42.9 kDa |
| Product Description | |
| Description | <p>A denatured recombinant protein with a N-Terminal His-tag and corresponding to the amino acids 1-340 of Human RASSF1</p> <p>Source: <i>E.coli</i></p> <p>Amino Acid Sequence: MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMSE PELIELRELA PAGRAGKGRT RLERANALRI ARGTACNPTR QLVPGRGHRF QPAGPATHTW CDLCGDFIWG VVRKGLQCAH CKFTCHYRCR ALVCLDCCGP RDLGWEPAVE RDTNVDEPVE WETPDLSQAE IEQKIKEYNA QINSNLFMSL NKDGSYTGFI KVQLKLVPRV SVPSSKKPPS LQDARRGPGR GTSVRRRTSF YLPKDAVKHL HVLSRTRARE VIEALLRKFL VVDDPRKFAL FERAERHGQV YLRKLLDDEQ PLRLRLLAGP SDKALSFLK ENDSGEVNW DAFSMPELHNF LRILQREEEE HLRQILQKYS YCRQKIQEAL HACPLG</p> |
| Gene ID | 11186 |
| Gene Symbol | RASSF1 |
| Species | Human |
| Product Application Details | |
| Applications | SDS-Page |
| Recommended Dilutions | SDS-Page |
| Application Notes | Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays. |

Images

SDS-Page: RASSF1 Protein [NBP2-23397]



15% SDS-PAGE (3ug)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-23397

| | |
|---------------|------------------------------------|
| NBP1-89411PEP | RASSF1 Recombinant Protein Antigen |
| AF748 | E-Cadherin Antibody [Unconjugated] |
| NBP2-03644 | RASSF1 Antibody (OTI2B11) |
| NB200-103 | p53 Antibody (PAb 240) - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-23397

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

