

Product Datasheet

Reg3A Antibody - BSA Free

NBP2-24763

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 3

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-24763

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-24763



NBP2-24763

Reg3A Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein G purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Rabbit Reg3A Antibody - BSA Free (NBP2-24763) is a polyclonal antibody validated for use in IHC and WB. Anti-Reg3A Antibody: Cited in 3 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5068
Gene Symbol	REG3A
Species	Human, Primate
Reactivity Notes	The amino acid sequence used as immunogen is 89% homologous in Monkey.
Immunogen	A portion of amino acids 1-50 of human REG3A was used as the immunogen for this antibody.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1-3 ug/ml, Immunohistochemistry, Immunohistochemistry-Paraffin 2.5-5 ug/ml




Images

Western Blot: Reg3A Antibody [NBP2-24763] - Human intestine lysate probed with REG3A antibody at 2 ug/ml. Goat anti-rabbit Ig HRP secondary antibody and PicoTect ECL substrate solution were used for this test.



Immunohistochemistry-Paraffin: Reg3A Antibody [NBP2-24763] - Human intestine stained with REG3A antibody (2.5 ug/ml), peroxidase-conjugate and DAB chromogen.



Publications

Jiang T, Huang Y, Liu H et al. Reduced miR-146a Promotes REG3A Expression and Macrophage Migration in Polymyositis and Dermatomyositis *Front Immunol* 2020-02-21 [PMID: 32153557] (IHC-P, Human)

Xu X, Fukui H, Ran Y et al. The Link between Type III Reg and STAT3-Associated Cytokines in Inflamed Colonic Tissues *Mediators Inflamm.* 2019-11-06 [PMID: 31780871] (WB, Human)

Obayashi N, Ohtsuka Y, Hosoi K et al. Comparison of Gene Expression Between Pediatric and Adult Gastric Mucosa with *Helicobacter pylori* Infection *Helicobacter* 2015-07-03 [PMID: 26140656] (IHC-P, Human)

Details:

Reg3A antibody was used for IHC-P staining on tissues from pediatric and adult human gastric mucosa with or without *Helicobacter pylori* infection.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-24763

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-24763

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

