

Product Datasheet

Nanog Antibody - BSA Free

NBP2-24941

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-24941

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-24941



NBP2-24941

Nanog Antibody - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein G purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Rabbit Nanog Antibody - BSA Free (NBP2-24941) is a polyclonal antibody validated for use in WB, Flow and ICC/IF. Anti-Nanog Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	79923
Gene Symbol	NANOG
Species	Human, Mouse
Marker	Embryonic Stem Cell Marker
Immunogen	A portion of amino acids 1-50 of human NANOG was used as the immunogen.


Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 0.5-2 ug/ml, Flow Cytometry reported in scientific literature (PMID 26283864), Immunocytochemistry/ Immunofluorescence 1:10

Images


Western Blot: Nanog Antibody [NBP2-24941] - Analysis of Nanog in (A) HeLa and (B) NIH 3T3 cell lysate using NBP2-24941 at 2 ug/ml.




Immunocytochemistry/Immunofluorescence: Nanog Antibody [NBP2-24941] - Nanog antibody was tested in Ntera2 cells with DyLight 488 (green). Nuclei and alpha-tubulin were counterstained with DAPI (blue) and DyLight 550 (red). Image objective 40x. An antibody dilution of 1:10 was used.




Western Blot: Nanog Antibody [NBP2-24941] - Protein expression of pluripotency markers in uterine cells isolated from control cows and from cows with adenomyosis. MM - molecular weight marker, C - cells obtained from control cows, ADENO - cells obtained from cows with adenomyosis. Image collected and cropped by CiteAb from the following publication (<https://rbej.biomedcentral.com/articles/10.1186/s12958-015-0106-0>), licensed under a CC-BY license.



Protein expression of NANOG a, OCT4 b & SOX2 c in bovine uterine tissues obtained from control cows & from cows with adenomyosis. Data were normalized against glyceraldehyde-3-phosphate dehydrogenase (GAPDH). Bars represent the mean \pm SEM. There were no statistical differences between uterine normal & adenomyotic tissues ($P > 0.05$), as determined by Student's t-test. Representative blots for NANOG, OCT4, SOX2 & GAPDH are shown below the graphs d. MM – molecular weight marker, C – tissues obtained from control cows, ADENO – tissues obtained from cows with adenomyosis Image collected & cropped by CiteAb from the following publication (<http://www.rbej.com/content/13/1/110>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Immunohistochemistry: Nanog Antibody [NBP2-24941] - Immunodetection of pluripotency markers in uterine tissues from control cows & from cows with adenomyosis. a, b – haematoxylin & eosin stained slides of control & adenomyotic uterus, respectively; c, d – NANOG immunodetection in control & adenomyotic tissue, respectively; e, f – immunolocalisation of OCT4 in normal & adenomyotic tissue, respectively; g, h – SOX2 immunodetection in control & adenomyotic tissue, respectively; i, j, k – no Ab, negative controls for NANOG, OCT4 & SOX2, respectively. Unspecific IgG controls (pictures not shown) served similar pictures as no Ab control. Arrows indicate the most intense histochemical reactions; dotted line indicate endometrial-myometrial border; e – endometrium, m – myometrium, es – endometrial stroma, eg – endometrial gland, a – adenomyotic lesion, v – vessel. Scale bars: 20 μ m Image collected & cropped by CiteAb from the following publication (<http://www.rbej.com/content/13/1/110>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Lupicka M, Socha B, Szczepanska A, Korzekwa A. Expression of pluripotency markers in the bovine uterus with adenomyosis. *Reprod Biol Endocrinol* 2015-01-01 [PMID: 26416515] (WB)

Dhamodaran K, Subramani M, Jeyabalan N et al. Characterization of ex vivo cultured limbal, conjunctival, and oral mucosal cells: A comparative study with implications in transplantation medicine. *Mol. Vis.* 2015-08-18 [PMID: 26283864] (FLOW, Human)

Details:

Nanog antibody was used for FLOW on limbal, conjunctival, and oral epithelial cells that were collected from human biopsy explant cultures established on the denuded human amniotic membrane with corneal lineage differentiation medium. The FLOW assay involved trypsinization of day 14 cultures, fixation using 4% PFA on ice for 10 minutes, permeabilization with 0.1% Triton X 100, primary detection with anti-rabbit Alexa 488 conjugated secondary antibody (Figure 1).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-24941

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-24941

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

