

Product Datasheet

Isocitrate Dehydrogenase 1/IDH1 Antibody - BSA Free NBP2-32150

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-32150

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-32150



NBP2-32150**Isocitrate Dehydrogenase 1/IDH1 Antibody - BSA Free**

Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)


Product Description	
Description	Novus Biologicals Rabbit Isocitrate Dehydrogenase 1/IDH1 Antibody - BSA Free (NBP2-32150) is a polyclonal antibody validated for use in WB and IP. Anti-Isocitrate Dehydrogenase 1/IDH1 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3417
Gene Symbol	IDH1
Species	Human, Mouse, Rat
Immunogen	The immunogen this antibody was made to, maps to a region between residue 364 to 414 of human Isocitrate Dehydrogenase 1 using the numbering given in entry NP_005887.2 (GeneID 3417).

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500 - 1:2500, Immunoprecipitation 2 - 10 ug/mg lysate




Images

Western Blot: Isocitrate Dehydrogenase 1/IDH1 Antibody [NBP2-32150]
 - Samples: Whole cell lysate (50 ug) from K562, HeLa, 293T, mouse TCMK-1, mouse NIH3T3, and rat C6 cells. Antibodies: Affinity purified rabbit anti-Isocitrate dehydrogenase antibody NBP2-32150 used for WB at 1 ug/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.



Immunoprecipitation: Isocitrate Dehydrogenase 1/IDH1 Antibody [NBP2-32150]
 - Samples: Whole cell lysate (1 mg for IP; 20% of IP loaded) from Jurkat cells. Antibodies: Affinity purified rabbit anti-Isocitrate dehydrogenase antibody NBP2-32150 used for IP at 6 ug/mg lysate. For blotting immunoprecipitated Isocitrate dehydrogenase, NBP2-32150 was used at 1 ug/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.



Publications

Malik RA, Zippel N, Fromel T et al. AMP-Activated Protein Kinase alpha2 in Neutrophils Regulates Vascular Repair via Hypoxia-Inducible Factor-1alpha and a Network of Proteins Affecting Metabolism and Apoptosis. *Circulation Research*. 2016-10-24 [PMID: 27777247] (WB, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-32150

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-32150

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

