

# Product Datasheet

## TLR5 Antibody (85B152.5) - Azide and BSA Free NBP2-33253

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-33253](http://www.novusbio.com/NBP2-33253)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-33253](http://www.novusbio.com/reviews/destination/NBP2-33253)



**NBP2-33253**

TLR5 Antibody (85B152.5) - Azide and BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	85B152.5
Preservative	No Preservative
Isotype	IgG2a
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	100 kDa

Product Description	
Description	Novus Biologicals Mouse TLR5 Antibody (85B152.5) - Azide and BSA Free (NBP1-97728) is a monoclonal antibody validated for use in WB and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	7100
Gene Symbol	TLR5
Species	Human, Mouse, Canine
Immunogen	This antibody was developed against KLH-conjugated synthetic peptide corresponding to a portion of human TLR5. It will cross-react with mouse TLR5 (NP_003259).


Product Application Details	
Applications	Western Blot, Flow Cytometry, Flow (Cell Surface)
Recommended Dilutions	Western Blot 1-3 ug/ml, Flow Cytometry, Flow (Cell Surface) reported in scientific literature (Wong et al (2007))
Application Notes	The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.

**Images**


Western Blot: TLR5 Antibody (85B152.5) - Azide Free [NBP2-33253] - Normal human bronchial epithelial total cell lysate. Image from verified customer review. Image from the standard format of this antibody.




Flow Cytometry: TLR5 Antibody (85B152.5) - Azide Free [NBP2-33253] - Intracellular flow analysis of TLR5 in  $10^6$  human lymphocytes using 0.5 ug of NB200-571. The shaded histogram represents cells without antibody, green represents isotype control (BD), and red represents TLR5 antibody. Image from the standard format of this antibody.




Flow Cytometry: TLR5 Antibody (85B152.5) - Azide Free [NBP2-33253] - Analysis using the FITC conjugate of NBP1-97728. Staining of TLR5 in  $10^6$  cells using 0.5 ugs of FITC-conjugated antibody: A) Intracellular flow testing of human Ramos cells and B) cell surface flow testing of mouse RAW cells. The shaded (or black) histogram represents cells without antibody, green represents isotype control antibody, and red represents TLR5 antibody.




Flow Cytometry: TLR5 Antibody (85B152.5) - Azide Free [NBP2-33253] - Analysis using the PE conjugate of NBP1-97728. Staining of TLR5 in  $10^6$  human lymphocytes using 0.5 ugs of PE-conjugated antibody. The shaded histogram represents cells without antibody, green represents isotype control antibody, and red represents TLR5 antibody.



Flow Cytometry: TLR5 Antibody (85B152.5) - Azide Free [NBP2-33253] - Analysis using the PE conjugate of NBP1-97728. Staining of TLR5 in  $10^6$  RAW cells using 0.5 ugs of PE-conjugated antibody. The black histogram represents cells without antibody, green represents isotype control antibody, and red represents TLR5 antibody.



Flow Cytometry: TLR5 Antibody (85B152.5) - Azide Free [NBP2-33253] - Intracellular flow analysis of TLR5 in Ramos cells using 1 ug of TLR5 antibody. Shaded histogram represents Ramos cells without antibody; green represents isotype control antibody; red represents TLR5 antibody. Image from the standard format of this antibody.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-33253**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-33253](http://www.novusbio.com/reviews/submit/NBP2-33253)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

