

Product Datasheet

Endophilin B1/Bif-1 Antibody (30A882.1.1) - Azide and BSA Free NBP2-33258

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-33258

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-33258



NBP2-33258


Endophilin B1/Bif-1 Antibody (30A882.1.1) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	30A882.1.1
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Protein G purified
Buffer	PBS
Product Description	
Description	Novus Biologicals Mouse Endophilin B1/Bif-1 Antibody (30A882.1.1) - Azide and BSA Free (NBP2-24733) is a monoclonal antibody validated for use in IHC, WB, ICC/IF, Simple Western and IP. Anti-Endophilin B1/Bif-1 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	51100
Gene Symbol	SH3GLB1
Species	Human, Mouse, Rat
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 25619828)
Immunogen	A full-length human Bif-1 recombinant protein was used as immunogen.
Product Application Details	
Applications	Western Blot, Simple Western, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1-2 ug/ml, Simple Western 1:200, Immunohistochemistry 5 ug/ml, Immunocytochemistry/ Immunofluorescence 1-2 ug/ml, Immunoprecipitation 1:10 -1:500, Immunohistochemistry-Paraffin 5 ug/ml
Application Notes	See Simple Western Antibody Database for Simple Western validation: antibody dilution of 1:200




Images

Western Blot: Endophilin B1/Bif-1 Antibody (30A882.1.1) - Azide Free [NBP2-33258] - Analysis for Bi-f-1/Endophilin-B1 using this antibody at 2 ug/ml dilution against 10 ug of Jurkat cell lysate.



Immunohistochemistry-Paraffin: Endophilin B1/Bif-1 Antibody (30A882.1.1) - Azide Free [NBP2-33258] - IHC analysis of formalin fixed paraffin-embedded (FFPE) human kidney using 5 ug/ml conc. of Endophilin B1 antibody (clone 30A882.1.1) on a Bond Rx autostainer (Leica Biosystems). The assay involved 20 minutes of heat induced antigen retrieval (HIER) using 10mM sodium citrate buffer (pH 6.0) and endogenous peroxidase quenching with peroxide block. The sections were incubated with primary antibody for 30 minutes and Bond Polymer Refine Detection (Leica Biosystems) with DAB was used for signal development followed by counterstaining with hematoxylin. Whole slide scanning and capturing of representative images was performed using Aperio AT2 (Leica Biosystems). Cytoplasmic staining of Endophilin B1 was observed. Staining was performed by Histowiz. Image using the standard format of this product.



Simple Western: Endophilin B1/Bif-1 Antibody (30A882.1.1) - Azide Free [NBP2-33258] - Simple Western lane view shows a specific band for Endophilin B1/Bif-1 in 0.05 mg/ml of Jurkat lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system.



Publications

Cho SG, Xiao X, Wang S et al Bif-1 Interacts with Prohibitin-2 to Regulate Mitochondrial Inner Membrane during Cell Stress and Apoptosis. *J Am Soc Nephrol*. 2019-07-01 [PMID: 31126972] (WB, IP, Human)

Details:

Citation using the Azide Free version of this antibody.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-33258

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-33258

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

