

Product Datasheet

CD40/TNFRSF5 Antibody - BSA Free NBP2-33957

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-33957

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-33957



NBP2-33957

CD40/TNFRSF5 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Rabbit CD40/TNFRSF5 Antibody - BSA Free (NBP2-33957) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	958
Gene Symbol	CD40
Species	Human
Immunogen	This antibody was developed against a recombinant protein corresponding to amino acids: EGWHCTSEACESCVLHRSCSPGFGVKQIATGVSDTICEPCPVGFFSNVSSAFE KCHPWTSCETKDLVVQQAGTNKTDVVCGPQD


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1.0 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunohistochemistry-Paraffin 1:50 - 1:200
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images



Western Blot: CD40/TNFRSF5 Antibody [NBP2-33957] - Lane 1: Marker [kDa] 250, 130, 100, 70, 55, 35, 25, 15, 10. Lane 2: U-2 OS




Immunohistochemistry-Paraffin: CD40/TNFRSF5 Antibody [NBP2-33957] - Staining of human tonsil shows strong membranous positivity in germinal center cells.




Immunohistochemistry-Paraffin: CD40/TNFRSF5 Antibody [NBP2-33957] - Analysis in human tonsil and cerebral cortex tissues. Corresponding CD40/TNFRSF5 RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: CD40/TNFRSF5 Antibody [NBP2-33957] - Staining of human cerebral cortex shows no positivity in neurons as expected.




Immunohistochemistry-Paraffin: CD40/TNFRSF5 Antibody [NBP2-33957] - Staining of human cerebral cortex, lymph node, spleen and tonsil using Anti-CD40/TNFRSF5 antibody NBP2-33957 (A) shows similar protein distribution across tissues to independent antibody NBP2-33956 (B).




Immunohistochemistry-Paraffin: CD40/TNFRSF5 Antibody [NBP2-33957] - Staining of human lymph node shows moderate to strong membranous positivity in germinal center cells.



Immunohistochemistry-Paraffin: CD40/TNFRSF5 Antibody [NBP2-33957] - Staining of human skeletal muscle shows no positivity in myocytes as expected.



Immunohistochemistry-Paraffin: CD40/TNFRSF5 Antibody [NBP2-33957] - Staining of human spleen shows moderate membranous positivity in cells in white pulp.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-33957

NBP2-33957PEP	CD40/TNFRSF5 Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-33957

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

