

# Product Datasheet

## enkurin Antibody - BSA Free

### NBP2-34037

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-34037](http://www.novusbio.com/NBP2-34037)

Updated 12/2/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-34037](http://www.novusbio.com/reviews/destination/NBP2-34037)



**NBP2-34037**

enkurin Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Rabbit enkurin Antibody - BSA Free (NBP2-34037) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	219670
Gene Symbol	ENKUR
Species	Human
Immunogen	This antibody was developed against a recombinant protein corresponding to amino acids: KKAQEDYDRYIQENLKKAAMKRLSDEEREAVLQGLKKNWEEVHKEFQSLSVFI DSIPKKIRKQRLEEEMKQLEHDIGIIEKHKIIYIAN


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:1000 - 1:2500, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

**Images**


Western Blot: enkurin Antibody [NBP2-34037] - Analysis in control (vector only transfected HEK293T lysate) and ENKUR over-expression lysate (Co-expressed with a C-terminal myc-DDK tag (3.1 kDa) in mammalian HEK293T cells).




Immunohistochemistry-Paraffin: enkurin Antibody [NBP2-34037] - Staining of human kidney.





Immunohistochemistry-Paraffin: enkurin Antibody [NBP2-34037] - Staining of human colon shows low expression as expected.




Immunohistochemistry-Paraffin: enkurin Antibody [NBP2-34037] - Staining of human fallopian tube shows high expression.




Immunohistochemistry-Paraffin: enkurin Antibody [NBP2-34037] - Staining in human fallopian tube and colon tissues using anti-ENKUR antibody. Corresponding ENKUR RNA-seq data are presented for the same tissues.



Immunohistochemistry-Paraffin: enkurin Antibody [NBP2-34037] - Staining of human bronchus, colon, fallopian tube and kidney using Anti-ENKUR antibody NBP2-34037 (A) shows similar protein distribution across tissues to independent antibody NBP1-89911 (B).



Immunohistochemistry-Paraffin: enkurin Antibody [NBP2-34037] - Staining of human bronchus.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-34037**

---

NBP2-34037PEP	enkurin Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-34037](http://www.novusbio.com/reviews/submit/NBP2-34037)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

