

# Product Datasheet

## KLHL33 Recombinant Protein Antigen NBP2-34165PEP

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-34165PEP](http://www.novusbio.com/NBP2-34165PEP)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-34165PEP](http://www.novusbio.com/reviews/destination/NBP2-34165PEP)



**NBP2-34165PEP****KLHL33 Recombinant Protein Antigen**

| <b>Product Information</b>     |                                                                                             |
|--------------------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>               | 0.1 ml                                                                                      |
| <b>Concentration</b>           | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>                 | Store at -20C. Avoid freeze-thaw cycles.                                                    |
| <b>Preservative</b>            | No Preservative                                                                             |
| <b>Purity</b>                  | >80% by SDS-PAGE and Coomassie blue staining                                                |
| <b>Buffer</b>                  | PBS and 1M Urea, pH 7.4.                                                                    |
| <b>Target Molecular Weight</b> | 27 kDa                                                                                      |

| <b>Product Description</b> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | <p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human KLHL33.</p> <p><b>Source:</b> <i>E. coli</i></p> <p><b>Amino Acid Sequence:</b><br/> LSPARCLALFPMAEAPGLERLWSKARHYLLTHLPAVALCPAFPSLPAACLAELL<br/> DSEDLHVQEEFEAFVAARCWLAANPETQES</p> <p><b>Fusion Tag:</b> N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G)</p> <p><b>This product is intended to be used as a blocking antigen for antibody competition assays. Any other use of this antigen is done at the risk of the user. The use of this product for commercial production is strictly prohibited. Please contact technical support if you have any questions.</b></p> |
| <b>Gene ID</b>             | 123103                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Gene Symbol</b>         | KLHL33                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Species</b>             | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |

| <b>Product Application Details</b> |                                            |
|------------------------------------|--------------------------------------------|
| <b>Applications</b>                | Antibody Competition                       |
| <b>Recommended Dilutions</b>       | Antibody Competition 10 - 100 molar excess |



**Application Notes**

This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP2-34165.

It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml.

For current lot information, including availability, please contact our technical support team click [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)

For further blocking peptide related information and a protocol, click [here](#).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-34165PEP**

---

|                  |                                    |
|------------------|------------------------------------|
| NBP3-17104PEP    | KLHL33 Recombinant Protein Antigen |
| NBP3-17104-100ul | KLHL33 Antibody                    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-34165PEP](http://www.novusbio.com/reviews/submit/NBP2-34165PEP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

