

Product Datasheet

E-Cadherin Antibody (CL1170) - BSA Free NBP2-34475

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-34475

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-34475



NBP2-34475

E-Cadherin Antibody (CL1170) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL1170
Preservative	0.02% Sodium Azide
Isotype	IgG2b
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Mouse E-Cadherin Antibody (CL1170) - BSA Free (NBP2-34475) is a monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	999
Gene Symbol	CDH1
Species	Human
Reactivity Notes	Rat (80%).
Immunogen	This antibody was developed using a recombinant protein derived from P12830, with the exact immunogen sequence remaining proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 1 ug/ml, Immunohistochemistry 1:200 - 1:500, Immunocytochemistry/ Immunofluorescence 2-10 ug/ml, Immunohistochemistry-Paraffin 1:200 - 1:500, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. For ICC/IF, PFA/Triton X-100 is recommended for fixation/permeabilization.

Images


Immunohistochemistry-Paraffin: E-Cadherin Antibody (CL1170) [NBP2-34475] - Staining in human duodenum and cerebral cortex tissues. Corresponding CDH1 RNA-seq data are presented for the same tissues.




Western Blot: E-Cadherin Antibody (CL1170) [NBP2-34475] - Analysis in MCF-7 cells transfected with control siRNA, target specific siRNA probe #1 and #2, using Anti-CDH1 antibody. Remaining relative intensity is presented. Loading control: Anti-GAPDH




Immunocytochemistry/Immunofluorescence: E-Cadherin Antibody (CL1170) [NBP2-34475] - Staining of BEWO cells using the Anti-CDH1 monoclonal antibody, showing specific staining in the plasma membrane and cell junctions in green. Microtubule- and nuclear probes are visualized in red and blue, respectively (where available).




Immunohistochemistry-Paraffin: E-Cadherin Antibody (CL1170) [NBP2-34475] - Staining of human duodenum shows strong membranous positivity in glandular cells.




Immunohistochemistry-Paraffin: E-Cadherin Antibody (CL1170) [NBP2-34475] - Staining of human prostate shows strong membranous positivity in glandular cells.



Immunohistochemistry-Paraffin: E-Cadherin Antibody (CL1170) [NBP2-34475] - Staining of human colorectal cancer shows moderate to strong membranous positivity in tumor cells.



Immunohistochemistry-Paraffin: E-Cadherin Antibody (CL1170) [NBP2-34475] - Staining of human cerebral cortex shows no positivity in neuronal cells as expected.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-34475

NBP2-34475PEP	E-Cadherin Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-34475

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

