

Product Datasheet

Siglec-3/CD33 Antibody (WM53) [PE/Cy7] NBP2-34511PECY7

Unit Size: 0.1 ml

Store at 4C in the dark. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-34511PECY7

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-34511PECY7



NBP2-34511PECY7

Siglec-3/CD33 Antibody (WM53) [PE/Cy7]

Product Information

Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark. Do not freeze.
Clonality	Monoclonal
Clone	WM53
Preservative	0.05% Sodium Azide
Isotype	IgG1 Kappa
Conjugate	PE/Cy7
Purity	Protein A or G purified
Buffer	PBS

Product Description


Host	Mouse
Gene ID	945
Gene Symbol	CD33
Species	Human
Immunogen	Human AML cells

Product Application Details

Applications	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, CyTOF-ready
Recommended Dilutions	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined. For optimal results using our Tandem dyes, please avoid prolonged exposure to light or extreme temperature fluctuations. These can lead to irreversible degradation or decoupling. When staining intracellular targets, specific attention to the fixation and permeabilization steps in your flow protocol may be required. Please contact our technical support team at technical@novusbio.com if you have any questions.

Images

Flow Cytometry: Siglec-3/CD33 Antibody (WM53) [PE/Cy7] [NBP2-34511PECY7] - Monocyte differentiation stages characterized by CD45, CD34, CD33 (NBP2-34511PECY7) and CD14 expression. HLA-DR and CCR2 expression were measured in each stage and expressed as the number of sites per cell. A) example of FACs analysis for monocyte differentiation (red: myelo/monoblasts, blue: promonocytes, green: monocytes), and the HLA-DR and CCR2 expression levels in each stage. Image collected and cropped by CiteAb from the following publication ([//pubmed.ncbi.nlm.nih.gov/27893741/](https://pubmed.ncbi.nlm.nih.gov/27893741/)) licensed under a CC-BY license.



Publications

Faivre V, Lukaszewicz A. C, et al. Downregulation of Blood Monocyte HLA-DR in ICU Patients Is Also Present in Bone Marrow Cells. PLoS One 2016-11-29 [PMID: 27893741] (FLOW, Human)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-34511PECY7

NBP2-32709PEP	Siglec-3/CD33 Recombinant Protein Antigen
DC140	CD14 [HRP]
1137-SL-050	Siglec-3/CD33 [Unconjugated]
AF3667	Myeloperoxidase/MPO Antibody [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-34511PECY7

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

