

Product Datasheet

NIFK Antibody (CL2237) - BSA Free NBP2-36750

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-36750

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-36750



NBP2-36750

NIFK Antibody (CL2237) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL2237
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Mouse NIFK Antibody (CL2237) - BSA Free (NBP2-36750) is a monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	84365
Gene Symbol	NIFK
Species	Human
Immunogen	This antibody was developed using a recombinant protein derived from Q9BYG3, with the exact immunogen sequence remaining proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1 ug/ml, Immunohistochemistry 1:2500 - 1:5000, Immunocytochemistry/ Immunofluorescence 2-10 ug/ml, Immunohistochemistry-Paraffin 1:2500 - 1:5000
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. For ICC/IF, PFA/Triton X-100 is recommended for fixation/permeabilization.

Images


Western Blot: NIFK Antibody (CL2237) [NBP2-36750] - Lane 1: Marker [kDa]. Lane 2: Human cell line U-251 MG




Immunocytochemistry/Immunofluorescence: NIFK Antibody (CL2237) [NBP2-36750] - Staining of U-2 OS cells using the Anti-MK167 monoclonal antibody, showing specific staining in the nucleoli in green. Microtubule- and nuclear probes are visualized in red and blue, respectively (where available).




Immunohistochemistry: NIFK Antibody (CL2237) [NBP2-36750] - Staining of human fallopian tube shows strong nucleolar immunoreactivity in glandular and connective tissue cells.




Immunohistochemistry: NIFK Antibody (CL2237) [NBP2-36750] - Staining of human skin shows strong nucleolar positivity in the epithelial cells.




Immunohistochemistry: NIFK Antibody (CL2237) [NBP2-36750] - Staining of human small intestine shows nucleolar immunoreactivity in glandular and lamina propria cells.




Immunohistochemistry: NIFK Antibody (CL2237) [NBP2-36750] - Staining of human pancreas shows strong nucleolar positivity in glandular cells.



Immunohistochemistry: NIFK Antibody (CL2237) [NBP2-36750] - Staining of human kidney shows strong nucleolar immunoreactivity in renal tubules.



Immunohistochemistry: NIFK Antibody (CL2237) [NBP2-36750] - Staining of human liver shows strong nucleolar positivity in glandular cells.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-36750

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-36750

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

