

Product Datasheet

PSP/BPIFA2 Antibody (3F8E4) - BSA Free NBP2-37277

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-37277

Updated 3/4/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-37277



NBP2-37277

PSP/BPIFA2 Antibody (3F8E4) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3F8E4
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	24.8 kDa

Product Description	
Host	Mouse
Gene ID	140683
Gene Symbol	BPIFA2
Species	Human
Immunogen	Purified recombinant fragment of mouse mSplunc2 (AA: 16-169) expressed in E. Coli.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Flow Cytometry, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Western Blot 1:500 - 1:2000, Flow Cytometry 1:200 - 1:400, ELISA 1:10000, Immunohistochemistry 1:200 - 1:1000, Immunohistochemistry-Paraffin 1:200 - 1:1000, CyTOF-ready
Application Notes	This antibody is Cytof ready.

Images


Western Blot: PSP/BPIFA2 Antibody (3F8E4) [NBP2-37277] - Western blot analysis using mSplunc2 mAb against HEK293 (1) and mSplunc2 - hlgGfc transfected HEK293 (2) cell lysate.




Immunohistochemistry: PSP/BPIFA2 Antibody (3F8E4) [NBP2-37277] - Immunohistochemical analysis of paraffin-embedded endometrial cancer tissues using mSplunc2 mouse mAb with DAB staining.




Flow Cytometry: PSP/BPIFA2 Antibody (3F8E4) [NBP2-37277] - Flow cytometric analysis of MCF-7 cells using mSplunc2 mouse mAb (green) and negative control (red).




Western Blot: PSP/BPIFA2 Antibody (3F8E4) [NBP2-37277] - Western blot analysis using mSplunc2 mAb against human mSplunc2 (AA:) recombinant protein. (Expected MW is 18.5 kDa)



Immunohistochemistry: PSP/BPIFA2 Antibody (3F8E4) [NBP2-37277] - Immunohistochemical analysis of paraffin-embedded ovarian cancer tissues using mSplunc2 mouse mAb with DAB staining.



ELISA: PSP/BPIFA2 Antibody (3F8E4) [NBP2-37277] - Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-37277

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP1-91704PEP	PSP/BPIFA2 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-37277

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

