

Product Datasheet

TLR6 Antibody (86B1153.2) - Azide and BSA Free NBP2-37704

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-37704

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-37704



NBP2-37704

TLR6 Antibody (86B1153.2) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	86B1153.2
Preservative	No Preservative
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Mouse TLR6 Antibody (86B1153.2) - Azide and BSA Free (NB100-56536) is a monoclonal antibody validated for use in IHC and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10333
Gene Symbol	TLR6
Species	Human
Immunogen	This antibody was developed against a synthetic peptide corresponding to amino acids 408-424 of human TLR6.

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Flow Cytometry, Flow (Cell Surface), Immunohistochemistry
Recommended Dilutions	Flow Cytometry, Immunohistochemistry 5 ug/ml, Immunohistochemistry-Paraffin 5 ug/ml, Flow (Cell Surface) 0.5-1.0 ug/million cells




Images

Immunohistochemistry: TLR6 Antibody (86B1153.2) - Azide Free [NBP2-37704] - Tonsil probed with TLR6 antibody at 5 ug/ml. Novus's human tissue TMA was used for this test. Staining of formalin-fixed tissues is enhanced by boiling tissue sections in 10 mM sodium citrate buffer, pH 6.0 for 10-20 min followed by cooling at RT for 20 m



Flow Cytometry: TLR6 Antibody (86B1153.2) - Azide Free [NBP2-37704] - Cell surface analysis of TLR6 antibody in 10^6 human monocytes using 0.5 ug of this antibody. Shaded histogram represents cells without antibody; green represents isotype control antibody ; red represents anti-TLR6 antibody. goat anti-mouse IgG PE conjugated secondary antibody was used. Image from the standard format of this antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-37704

NBP2-26219-2ug	MALP-2, TLR6 and TLR2 ligand
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-37704

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

