

Product Datasheet

GAPDH Antibody (GA1R) - Azide and BSA Free NBP2-37828

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-37828

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-37828




NBP2-37828

GAPDH Antibody (GA1R) - Azide and BSA Free


Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	GA1R
Preservative	No Preservative
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2)
Target Molecular Weight	36 kDa
Product Description	
Description	Novus Biologicals Mouse GAPDH Antibody (GA1R) - Azide and BSA Free (NBP2-37828) is a monoclonal antibody validated for use in IHC, WB, ELISA, Flow, ICC/IF and Simple Western. Anti-GAPDH Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	2597
Gene Symbol	GAPDH
Species	Human, Mouse, Rat, Porcine, Bacteria, Chicken, Hamster, Insect, Rabbit, Yeast
Reactivity Notes	Porcine reactivity reported in scientific literature (PMID: 17868880). Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Additional Mouse on Mouse blocking steps may be required for IHC and ICC experiments. Please contact Technical Support for more information.
Specificity/Sensitivity	This GAPDH antibody detects GAPDH from BL-21 bacteria, Sf9 insect, Saccharomyces cerevisiae (yeast), human, mouse, rat, rabbit, hamster, and chicken samples. This product has been successfully used in Western blot, ELISA, and immunofluorescence.
Immunogen	This GAPDH protein was developed against recombinant GAPDH
Product Application Details	
Applications	Western Blot, Simple Western, Immunohistochemistry-Paraffin, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Western Blot 1:1000 - 1:10000, Simple Western, Flow Cytometry 2 ug, ELISA 1:100 - 1:2000, Immunohistochemistry 1:50 - 1:500, Immunocytochemistry/ Immunofluorescence 1:20 - 1:200, Immunohistochemistry-Paraffin 1:50 - 1:500, CyTOF-ready
Application Notes	See Simple Western Antibody Database for Simple Western validation: separated by Size

Images


Western Blot: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis of BL-21 bacteria (1), Sf9 insect (2), and *Saccharomyces cerevisiae* (3).
Theoretical molecular weight: 36 kDa.




Immunocytochemistry/Immunofluorescence: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis of GAPDH (green) in HeLa. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 10 minutes at room temperature and blocked with 5% normal goat serum for 15 minutes at room temperature. Cells were probed with a GAPDH monoclonal antibody at a dilution of 1:100 for at least 1 hour at room temperature, washed with PBS, and incubated with DyLight 488 goat anti-mouse IgG secondary antibody at a dilution of 1:400 for 30 minutes at room temperature. Nuclei (blue) were stained with Hoechst 33342 dye.




Immunohistochemistry-Paraffin: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis showing staining in the nucleus and cytoplasm of paraffin-embedded human breast carcinoma (right) compared with a negative control without primary antibody (left).




Flow Cytometry: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis of GAPDH in HeLa cells (green) compared to an isotype control (blue).




Western Blot: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis of 5ug/lane tissue lysates from human (1), mouse (2), rat (3), rabbit (4), chicken (5), and hamster (6). Theoretical molecular weight: 36 kDa.




Western Blot: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis of 20ug of THP-1 whole cell lysates per well onto an SDS-PAGE gel. GAPDH theoretical molecular weight: 36 kDa.




Immunocytochemistry/Immunofluorescence: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis of GAPDH (green) showing staining in the in the cytoplasm of MCF-7 cells (right) compared to a negative control without primary antibody (left).




Immunocytochemistry/Immunofluorescence: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis of GAPDH (green) showing staining in the in the cytoplasm of Hela cells (right) compared to a negative control without primary antibody (left).




Immunocytochemistry/Immunofluorescence: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis of GAPDH (green) showing staining in the in the cytoplasm of C2C12 cells (right) compared to a negative control without primary antibody (left).




Immunohistochemistry-Paraffin: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis showing staining in the nucleus and cytoplasm of mouse colon tissue (right) compared with a negative control without primary antibody (left).




Immunohistochemistry-Paraffin: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis showing staining in the nucleus and cytoplasm of human kidney tissue (right) compared with a negative control without primary antibody (left).



Flow Cytometry: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis of GAPDH in NIH-3T3 cells (green) compared to an isotype control (blue).



Flow Cytometry: GAPDH Antibody (GA1R) [NBP2-37828] - Analysis of GAPDH in MCF-7 cells (green) compared to an isotype control (blue).



Publications

Shi Z, Ren M, Rockey DC Myocardin and myocardin-related transcription factor-A synergistically mediate actin cytoskeletal dependent inhibition of liver fibrogenesis *Am. J. Physiol. Gastrointest. Liver Physiol.* 2020-01-13 [PMID: 31928221]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-37828

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-37828

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

