

# Product Datasheet

## TRACP/PAP/ACP5 Antibody [PE]

### NBP2-41178PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-41178PE](http://www.novusbio.com/NBP2-41178PE)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-41178PE](http://www.novusbio.com/reviews/destination/NBP2-41178PE)



**NBP2-41178PE**

TRACP/PAP/ACP5 Antibody [PE]

| Product Information         |                                                                                                                                                                                                                  |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                           |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                      |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                         |
| Clonality                   | Polyclonal                                                                                                                                                                                                       |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                               |
| Isotype                     | IgG                                                                                                                                                                                                              |
| Conjugate                   | PE                                                                                                                                                                                                               |
| Purity                      | Peptide affinity purified                                                                                                                                                                                        |
| Buffer                      | PBS                                                                                                                                                                                                              |
| Product Description         |                                                                                                                                                                                                                  |
| Host                        | Rabbit                                                                                                                                                                                                           |
| Gene ID                     | 54                                                                                                                                                                                                               |
| Gene Symbol                 | ACP5                                                                                                                                                                                                             |
| Species                     | Human, Mouse                                                                                                                                                                                                     |
| Reactivity Notes            | Immunogen displays the following percentage of sequence identity for non-tested species: Rat (89%), Rabbit (88%)                                                                                                 |
| Immunogen                   | Antibody was raised against an 18 amino acid synthetic peptide near the carboxy terminus of human TRAP. The immunogen is located within the last 50 amino acids of TRAP. Amino Acid Sequence: EASGKSLFKTRLPRRARP |
| Product Application Details |                                                                                                                                                                                                                  |
| Applications                | Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence                                                                                                                                                     |
| Recommended Dilutions       | Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence                                                                                                                                                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                           |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-41178PE**

---

|               |                                            |
|---------------|--------------------------------------------|
| NBP2-24983    | Rabbit IgG Isotype Control [PE]            |
| NBP2-38833PEP | TRACP/PAP/ACP5 Recombinant Protein Antigen |
| 210-TA-005    | TNF-alpha [Unconjugated]                   |
| 3948-AP-010   | TRACP/PAP/ACP5 [Unconjugated]              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-41178PE](http://www.novusbio.com/reviews/submit/NBP2-41178PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

