

Product Datasheet

IL-36Ra/IL-1F5 Antibody - BSA Free NBP2-41241

Unit Size: 0.1 mg

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-41241

Updated 3/4/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-41241



NBP2-41241

IL-36Ra/IL-1F5 Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	24 kDa


Product Description	
Host	Rabbit
Gene ID	26525
Gene Symbol	IL36RN
Species	Human, Mouse
Specificity/Sensitivity	At least two isoforms of IL-36RN are known to exist
Immunogen	Antibody was raised against a 19 amino acid peptide near the carboxy terminus of human IL-36RN. The immunogen is located within the last 50 amino acids of IL-36RN. Amino Acid Sequence: PEADQPVRTQLPENGAWN

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1 ug/ml, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 20 ug/ml, Immunohistochemistry-Paraffin 1:10-1:500




Images


Western Blot: IL-36Ra/IL-1F5 Antibody [NBP2-41241] - Analysis of IL-36RN in human spleen tissue lysate with IL-36RN antibody at 1 ug/ml in (A) the absence and (B) the presence of blocking peptide.




Immunohistochemistry: IL-36Ra/IL-1F5 Antibody [NBP2-41241] - Immunohistochemistry of IL-36RN in rat small intestine tissue with IL-36RN antibody at 5 ug/mL.



Flow Cytometry: IL-36Ra/IL-1F5 Antibody [NBP2-41241] - Het1a normal esophageal cell line. Basal expression of IL-36a receptor on the surface of normal esophageal cells. Flow cytometry image submitted by a verified customer review. Image from the PE version of this antibody.



Immunocytochemistry/ Immunofluorescence: IL-36Ra/IL-1F5 Antibody - BSA Free [NBP2-41241] - Immunofluorescence of IL-36Ra/IL-1F5 in rat small intestine tissue with IL-36Ra/IL-1F5 antibody at 20 ug/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-41241

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-34872-100ug	Recombinant Human IL-36Ra/IL-1F5 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-41241

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

