

Product Datasheet

TAF9 Antibody - BSA Free NBP2-41266

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-41266

Updated 3/4/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-41266



NBP2-41266

TAF9 Antibody - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	28 kDa

Product Description	
Host	Rabbit
Gene ID	6880
Gene Symbol	TAF9
Species	Human, Mouse, Rat
Specificity/Sensitivity	Multiple isoforms of TAF9 are known to exist. TAF9 is predicted to not cross-react with other TAF family members.
Immunogen	Antibody was raised against a 17 amino acid peptide near the center of human TAF9. The immunogen is located within amino acids 110 - 160 of TAF9. Amino Acid Squence: LTAPNYRLKSLQKKAS


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1-2 ug/ml, ELISA, Immunohistochemistry 5 ug/ml, Immunocytochemistry/ Immunofluorescence 20 ug/ml, Immunohistochemistry-Paraffin 5 ug/ml

Images


Immunohistochemistry: TAF9 Antibody [NBP2-41266] - TAF9 in mouse liver tissue with TAF9 antibody at 5 ug/mL.




Immunohistochemistry: TAF9 Antibody - BSA Free [NBP2-41266] - Immunohistochemistry of TAF9 in mouse liver tissue with TAF9 antibody at 5 ug/mL.




Immunocytochemistry/ Immunofluorescence: TAF9 Antibody - BSA Free [NBP2-41266] - Immunofluorescence of TAF9 in mouse liver tissue with TAF9 antibody at 20 u/mL.



Immunohistochemistry: TAF9 Antibody - BSA Free [NBP2-41266] - Immunohistochemistry of TAF9 in mouse liver tissue with TAF9 antibody at 5 u/mL.



Immunocytochemistry/ Immunofluorescence: TAF9 Antibody - BSA Free [NBP2-41266] - Immunofluorescence of TAF9 in mouse liver tissue with TAF9 antibody at 20 u/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-41266

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-50978-0.1mg	Recombinant Human TAF9 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-41266

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

