

Product Datasheet

alpha-Synuclein Antibody - BSA Free NBP2-41288

Unit Size: 0.1 mg

Antibody can be stored at 4C up to one year. Antibodies should not be exposed to prolonged high temperatures.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-41288

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-41288



NBP2-41288

alpha-Synuclein Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Antibody can be stored at 4C up to one year. Antibodies should not be exposed to prolonged high temperatures.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	16 kDa


Product Description	
Description	Novus Biologicals Rabbit alpha-Synuclein Antibody - BSA Free (NBP2-41288) is a polyclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	6622
Gene Symbol	SNCA
Species	Human, Mouse, Rat
Immunogen	This alpha-Synuclein Antibody was developed against 19 amino acid peptide near the carboxy terminus of human alpha-Synuclein. The immunogen is located within the last 50 amino acids of alpha-Synuclein. Amino Acid Sequence: DMPVDPDNEAYEMPSEEGY

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1-2 ug/ml, ELISA, Immunohistochemistry 5 ug/ml, Immunocytochemistry/ Immunofluorescence 20 ug/ml, Immunohistochemistry-Paraffin 5 ug/ml




Images


Western Blot: alpha-Synuclein Antibody [NBP2-41288] - Analysis of 1 ng of recombinant human SNCA with SNCA antibody at 1 ug/ml.




Western Blot: alpha-Synuclein Antibody [NBP2-41288] - Western blot analysis of SNCA in mouse cerebellum tissue lysate with SNCA antibody at 1 ug/ml.



Immunohistochemistry: alpha-Synuclein Antibody - BSA Free [NBP2-41288] - Immunohistochemistry of alpha-Synuclein in rat brain tissue with alpha-Synuclein antibody at 5 ug/ml.



Immunocytochemistry/ Immunofluorescence: alpha-Synuclein Antibody - BSA Free [NBP2-41288] - Immunofluorescence of alpha-Synuclein in rat brain tissue with alpha-Synuclein antibody at 20 ug/ml.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-41288

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-54789-100ug	Recombinant Human alpha-Synuclein Active, Pre-formed Fibrils, (Type 1) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-41288

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

