

Product Datasheet

APC-PCI Complex Antibody (M36) - BSA Free NBP2-41380

Unit Size: 100 ug

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-41380

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-41380



NBP2-41380

APC-PCI Complex Antibody (M36) - BSA Free


Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	M36
Preservative	0.09% Sodium Azide
Isotype	IgG1 Kappa
Purity	Protein A or G purified
Buffer	0.01 M phosphate buffer, pH 7.4, with 0.14 M NaCl

Product Description	
Description	Novus Biologicals Mouse APC-PCI Complex Antibody (M36) - BSA Free (NBP2-41380) is a monoclonal antibody validated for use in ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Species	Human
Specificity/Sensitivity	Antibody is specific for APC:PCI complex and cleaved PCI but has little reactivity towards native PCI. The dissociation constants (as determined by standard surface plasmon resonance technique) are: $K_d = 4 \times 10^{-10}$ M for PCI in complex with APC and 2×10^{-10} M for cleaved PCI. The value of K_d for native PCI is too low to measure, i.e. $> 10^{-5}$ M.
Immunogen	Mixture of PCI in complex with PSA (approx. 80%) and PCI cleaved from such a complex (approx. 20%).

Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA 1:100-1:2000
Application Notes	ELISA: NBP2-41380 shows excellent reactivity towards PCI in cleaved and complex form in ELISA. NBP2-41380 is a part of the APC-PCI Matched Reagent Set, and can be used in ELISA as a capture antibody in combination with SA001RA as the detection antibody. (1-5) WB: NBP2-41380 shows low reactivity against denatured antigen in SDS-PAGE.

Images

ELISA: APC-PCI Complex Antibody (M36) - BSA Free [NBP2-41380] - The figure illustrates the calibration curve of a sandwich assay with NBP2-41380 as the capture antibody and biotinylated detection antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-41380

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-41380

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

