

Product Datasheet

HABP1/C1QBP/GC1q R Antibody (1G7) - BSA Free NBP2-42691

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-42691

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-42691



NBP2-42691

HABP1/C1QBP/GC1q R Antibody (1G7) - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1G7
Preservative	0.02% Sodium Azide
Isotype	IgG1 Kappa
Purity	Protein A purified
Buffer	PBS (pH 7.4), 10% Glycerol


Product Description	
Description	Novus Biologicals Mouse HABP1/C1QBP/GC1q R Antibody (1G7) - BSA Free (NBP2-42691) is a monoclonal antibody validated for use in WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	708
Gene Symbol	C1QBP
Species	Human
Immunogen	Recombinant human HABP1/C1QBP/GC1q R (74-282aa) purified from E.coli.

Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:2000, ELISA, Immunocytochemistry/ Immunofluorescence 1:100




Images

Western Blot: HABP1/C1QBP/GC1q R Antibody (1G7) [NBP2-42691] - The recombinant protein(50ng) and cell lysates(40ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human C1QBP antibody (1:2000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1 : Recombinant protein Lane 2 : HeLa cell lysate Lane 3 : Jurkat cell lysate Lane 4 : MCF7 cell lysate Lane 5 : Ramos cell lysate Lane 6 : 293T cell lysate Lane 7 : A549 cell lysate



Immunocytochemistry/Immunofluorescence: HABP1/C1QBP/GC1q R Antibody (1G7) [NBP2-42691] - Analysis of C1QBP in HeLa cells line, stained with DAPI (Blue) for nucleus staining and monoclonal anti-human C1QBP antibody (1:100) with goat anti-mouse IgG-Alexa fluor 488 conjugate (Green).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-42691

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-42691

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

