

Product Datasheet

SNRPE Antibody NBP2-43792

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-43792

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-43792




NBP2-43792

SNRPE Antibody


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 1% BSA, 20% Glycerol
Target Molecular Weight	11 kDa
Product Description	
Description	Novus Biologicals Rabbit SNRPE Antibody (NBP2-43792) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-SNRPE Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	6635
Gene Symbol	SNRPE
Species	Human, Mouse
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the center region of human SNRPE. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry Assay dependent, Immunocytochemistry/ Immunofluorescence Reported in scientific literature (PMID:34096600), Immunohistochemistry-Paraffin Assay dependent

Images


Western Blot: SNRPE Antibody [NBP2-43792] - Various whole cell extracts (30 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with SNRPE antibody diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody



SNRPE antibody detects SNRPE protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded mouse brain. SNRPE stained by SNRPE antibody (NBP2-43792) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



SNRPE antibody detects SNRPE protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded mouse brain. SNRPE stained by SNRPE antibody (NBP2-43792) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Publications

Esser LM, Schmitz K, Hillebrand F et al. Phosphorylation of pICln by the autophagy activating kinase ULK1 regulates snRNP biogenesis and splice activity of the cell Computational and structural biotechnology journal 2023-03-16 [PMID: 36968021] (WB, ICC/IF, Human)

Schmitz K, Cox J, Esser LM et al. An essential role of the autophagy activating kinase ULK1 in snRNP biogenesis Nucleic acids research 2021-06-21 [PMID: 34096600] (ICC/IF, WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-43792

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-43792

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

