

Product Datasheet

CXCR2/IL-8RB Antibody (547) - Azide and BSA Free NBP2-43810

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-43810

Updated 3/4/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-43810



NBP2-43810


CXCR2/IL-8RB Antibody (547) - Azide and BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	547
Preservative	No Preservative
Isotype	IgG2a
Purity	Protein A purified
Buffer	PBS, 20% Glycerol
Target Molecular Weight	41 kDa
Product Description	
Host	Mouse
Gene ID	3579
Gene Symbol	CXCR2
Species	Human
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of human CXCR2/IL-8RB. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:1000-1:10000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000




Images

Immunocytochemistry/Immunofluorescence: CXCR2/IL-8RB Antibody (547) [NBP2-43810] - CXCR2 antibody [547] detects CXCR2 protein by immunofluorescent analysis. Sample: Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: CXCR2 stained by CXCR2 antibody [547] diluted at 1:2000. Blue: Hoechst 33342 staining. Scale bar= 10um.



Non-transfected (-) and transfected (+) Boiled and unboiled 293T whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with CXCR2 antibody [GT547] (NBP2-43810) diluted at 1:5000. The HRP-conjugated anti-mouse IgG antibody was used to detect the primary antibody.



Publications

H Jayatilaka, P Tyle, JJ Chen, M Kwak, J Ju, HJ Kim, JSH Lee, PH Wu, DM Gilkes, R Fan, D Wirtz Synergistic IL-6 and IL-8 paracrine signalling pathway infers a strategy to inhibit tumour cell migration Nat Commun, 2017-05-26;8 (0):15584. 2017-05-26 [PMID: 28548090] (Western Blot, Mouse)

Jayatilaka H, Tyle P, Chen JJ et al. Synergistic IL-6 and IL-8 paracrine signalling pathway infers a strategy to inhibit tumour cell migration. Nat Commun. 2017 May 26 [PMID: 28548090] (WB, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-43810

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
H00003579-P01-10ug	Recombinant Human CXCR2/IL-8RB GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-43810

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

