

Product Datasheet

TGF-beta Antibody (1D11.16.8) - BSA Free NBP2-45137-0.1mg

Unit Size: 0.1 mg

Store at 4C short term. Store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 4

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-45137

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-45137



NBP2-45137-0.1mg

TGF-beta Antibody (1D11.16.8) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1D11.16.8
Preservative	0.02% Sodium Azide
Isotype	IgG1 Kappa
Purity	Protein A or G purified
Buffer	PBS
Target Molecular Weight	13 kDa


Product Description	
Description	Novus Biologicals Mouse TGF-beta Antibody (1D11.16.8) - BSA Free (NBP2-45137) is a monoclonal antibody validated for use in IHC, ELISA and ICC/IF. Anti-TGF-beta Antibody: Cited in 4 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	7040
Gene Symbol	TGFB1
Species	Human, Mouse, Rat, Bovine, Canine, Goat, Hamster, Monkey
Reactivity Notes	Goat reactivity reported from a verified customer review.
Specificity/Sensitivity	This MAb recognizes TGF beta 1, 2 and 3. Three TGF betas have been identified in mammals. TGFbeta1, TGFbeta2 and TGFbeta3 are each synthesized as precursor proteins that are very similar in that each is cleaved to yield a 112 amino acid polypeptide that remains associated with the latent portion of the molecules. Biologically active TGFbeta requires dimerization of the monomers (usually homodimers) and release of the latent peptide portion. Overall, the mature region of the TGFbeta3 protein has approximately 80% identity to the mature region of both TGFbeta1 and TGFbeta2. However, the NH2 terminals or precursor regions of their molecules share only 27% sequence identity. TGFbeta's inhibit the growth of epithelial cells and stimulate the growth of mesenchymal cells.
Immunogen	Bovine bone-derived TGF-beta 1 and TGF-beta 2

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Sandwich ELISA, ELISA Capture (Matched Antibody Pair)
Recommended Dilutions	Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 8-25 ug/ml, Immunohistochemistry-Paraffin 0.5 - 1.0ug/ml, Sandwich ELISA, ELISA Capture (Matched Antibody Pair)
Application Notes	Hu-chromosome location: 19q13.1 (beta1); 1q41 (beta2); 14q24 (beta3) For Cytokine Neutralization and Functional Assay, order Azide and BSA Free formulation (NBP2-47736)




Images


Immunocytochemistry/Immunofluorescence: TGF-beta Antibody (1D11.16.8) - BSA Free [NBP2-45137] - U2OS cells were fixed in 4% paraformaldehyde for 10 minutes and permeabilized in 0.05% Triton X-100 in PBS for 5 minutes. The cells were incubated with TGF-beta Antibody [1D11.16.8] (NBP2-45137) at 1ug/ml overnight at 4C and detected with an anti-mouse DyLight 488 (Green) at a 1:1000 dilution for 60 minutes. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 40X objective.



Immunohistochemistry: TGF-beta Antibody (1D11.16.8) [NBP2-45137] - TGF-beta was detected in immersion fixed paraffin-embedded human skin using 25 ug/mL Mouse Anti-TGF-beta Monoclonal Antibody (Clone 1D11.16.8) overnight at 4 degrees C. Tissues was stained using the Anti-Mouse HRP-DAB Cell and Tissue Staining Kit (brown; Catalog # CTS002) and counterstained with Haematoxylin (blue).



The sections of histological goat bone fragments were incubated with the primary antibody TGF-beta NBP2-45137 overnight. Subsequently, the secondary antibody NB7539 1:200; was applied and incubated for 1 hour. Image from a verified customer review.



Publications

Chen W, Pilling D, Gomer RH The mRNA binding protein DDX3 mediates TGF- β 1 upregulation of translation and promotes pulmonary fibrosis JCI insight 2023-02-23 [PMID: 36821384] (IHC-Fr, Mouse)

limuro Y, Yada A, Okada T et al. Cytoglobin-expressing cells in the splenic cords contribute to splenic fibrosis in cirrhotic patients Histol. Histopathol. 2020-09-18 [PMID: 32945524] (IF/IHC, IHC-P, Human)

Tashita C, Hoshi M, Hirata A et Al. Kynurenine plays an immunosuppressive role in 2,4,6-trinitrobenzene sulfate-induced colitis in mice World J. Gastroenterol. 2020-03-07 [PMID: 32206003] (IF/IHC, Mouse)

Scimeca M, Urbano N, Bonfiglio R et al. Breast osteoblast-like cells: a new biomarker for the management of breast cancer. Br. J. Cancer. 2018-10-17 [PMID: 30327566] (IHC-P, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-45137-0.1mg

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
7754-BH-005/CF	TGF-beta 1 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-45137

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

