

Product Datasheet

MBP Antibody (CL2829) - BSA Free NBP2-46632

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-46632

Updated 3/4/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-46632



NBP2-46632

MBP Antibody (CL2829) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL2829
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Host	Mouse
Gene ID	4155
Gene Symbol	MBP
Species	Human, Mouse, Rat
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions
Immunogen	This antibody was developed using a recombinant protein derived from P02686, with the exact immunogen sequence remaining proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry
Recommended Dilutions	Western Blot 1 ug/ml, ELISA Reported in scientific literature (PMID 29310668), Immunohistochemistry 1:5000 - 1:10000, Immunohistochemistry-Paraffin 1:5000 - 1:10000
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Immunohistochemistry-Paraffin: MBP Antibody (CL2829) [NBP2-46632] - Analysis in human cerebral cortex and liver tissues using NBP2-46632 antibody. Corresponding MBP RNA-seq data are presented for the same tissues.




Western Blot: MBP Antibody (2829) [NBP2-46632] - Lane 1: Marker [kDa] 250, 130, 100, 70, 55, 35, 25, 15, 10. Lane 2: Mouse Cerebral Cortex tissue




Immunohistochemistry-Paraffin: MBP Antibody (CL2829) [NBP2-46632] - Staining of human liver shows no positivity in hepatocytes as expected.




Western Blot: MBP Antibody (2829) [NBP2-46632] - Lane 1: Marker [kDa] 250, 130, 100, 70, 55, 35, 25, 15, 10. Lane 2: Human Cerebral Cortex tissue




Immunohistochemistry: MBP Antibody (2829) [NBP2-46632] - Staining of rat cerebral cortex shows strong immunoreactivity in myelinated fibers.




Immunohistochemistry: MBP Antibody (2829) [NBP2-46632] - Staining of mouse hippocampal formation shows strong immunoreactivity in myelinated fibers in the dentate gyrus.




Immunohistochemistry: MBP Antibody (2829) [NBP2-46632] - Staining of rat basal forebrain shows strong positivity in myelinated fibers.




Immunohistochemistry: MBP Antibody (2829) [NBP2-46632] - Staining of mouse cerebellum shows strong immunoreactivity in myelinated fibers.




Immunohistochemistry-Paraffin: MBP Antibody (2829) [NBP2-46632] - Staining of human cerebral cortex shows strong positivity in myelinated fibers.



Immunohistochemistry-Paraffin: MBP Antibody (CL2829) [NBP2-46632] - Staining of human cerebellum shows strong positivity in myelinated neural fibers.



Immunohistochemistry-Paraffin: MBP Antibody (CL2829) [NBP2-46632] - Staining of human pancreas shows no positivity in exocrine glandular cells as expected.



Publications

Brescia AC, Simonds MM, McCahan SM et al. Prior to extension, Transcriptomes of fibroblast-like Synoviocytes from extended and Polyarticular juvenile idiopathic arthritis are indistinguishable. *Pediatr Rheumatol Online J* 2018-01-08 [PMID: 29310668] (ELISA, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-46632

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP2-23253	Recombinant Human MBP His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-46632

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

