

Product Datasheet

Tyrosine Hydroxylase Antibody (CL3049) - BSA Free NBP2-46647

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-46647

Updated 3/4/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-46647



NBP2-46647

Tyrosine Hydroxylase Antibody (CL3049) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL3049
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol
Target Molecular Weight	60 kDa

Product Description	
Host	Mouse
Gene ID	7054
Gene Symbol	TH
Species	Human, Mouse, Rat
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions
Immunogen	This antibody was developed using a recombinant protein derived from P07101, with the exact immunogen sequence remaining proprietary.


Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 1:500 - 1:1000, Immunohistochemistry-Paraffin 1:500 - 1:1000
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Immunohistochemistry: Tyrosine Hydroxylase Antibody (CL3049) [NBP2-46647] - Staining of mouse brain shows strong positivity in neurons in the hypothalamus.




Immunohistochemistry: Tyrosine Hydroxylase Antibody (3049) [NBP2-46647] - Staining of human caudate nucleus shows immunoreactivity in monoaminergic neural fibers.




Immunohistochemistry: Tyrosine Hydroxylase Antibody (3049) [NBP2-46647] - Staining of human prostate shows strong immunoreactivity in sympathetic fibers.




Immunohistochemistry: Tyrosine Hydroxylase Antibody (3049) [NBP2-46647] - Analysis of human kidney shows absence of immunoreactivity (negative control).



Immunohistochemistry: Tyrosine Hydroxylase Antibody (3049) [NBP2-46647] - Immunohistochemical staining of rat locus coeruleus using Tyrosine Hydroxylase antibody (clone 3049) with immunofluorescence based detection showing a very strong cytoplasmic immunoreactivity in monoaminergic neurons.



Immunohistochemistry: Tyrosine Hydroxylase Antibody (CL3049) [NBP2-46647] - Staining of rat locus coeruleus shows strong positivity in noradrenergic neurons.



Publications

Kim S, Choi J, Kim S et al. Inhibition of alpha-synuclein aggregation by MT101-5 is neuroprotective in mouse models of Parkinson's disease *Biomedicine & Pharmacotherapy* 2022-10-01 [PMID: 36058149] (IHC-Fr, Mouse)

Details:
Brain Sections, dilution used 1:1000



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-46647

NBP2-46647PEP	Tyrosine Hydroxylase Recombinant Protein Antigen
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-46647

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

