

Product Datasheet

LMP7/PSMB8 Antibody - BSA Free NBP2-47395

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-47395

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-47395



NBP2-47395

LMP7/PSMB8 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Rabbit LMP7/PSMB8 Antibody - BSA Free (NBP2-47395) is a polyclonal antibody validated for use in IHC and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5696
Gene Symbol	PSMB8
Species	Human
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: Mouse (88%), Rat (84%)
Immunogen	This antibody was developed against a recombinant protein corresponding to amino acids: RPNLSPEEAYDLGRRAIAYATHRDS


Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 1:500 - 1:1000, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:500 - 1:1000
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. Immunocytochemistry/Immunofluorescence Fixation Permeabilization: Use PFA/Triton X-100.

Images


Immunocytochemistry/Immunofluorescence: LMP7/PSMB8 Antibody [NBP2-47395] - Immunofluorescent staining of human cell line HeLa shows localization to cytokinetic bridge & vesicles.




Immunohistochemistry-Paraffin: LMP7/PSMB8 Antibody [NBP2-47395] - Staining of human colon.




Immunohistochemistry-Paraffin: LMP7/PSMB8 Antibody [NBP2-47395] - Staining of human rectum shows moderate cytoplasmic and nuclear positivity in glandular cells.




Immunohistochemistry-Paraffin: LMP7/PSMB8 Antibody [NBP2-47395] - Staining of human lymph node shows high expression.




Immunohistochemistry-Paraffin: LMP7/PSMB8 Antibody [NBP2-47395] - Staining of human pancreas shows low expression as expected.




Immunohistochemistry-Paraffin: LMP7/PSMB8 Antibody [NBP2-47395] - Staining in human lymph node and pancreas tissues using anti-PSMB8 antibody. Corresponding PSMB8 RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: LMP7/PSMB8 Antibody [NBP2-47395] - Staining of human colon, kidney, liver and lymph node using Anti-PSMB8 antibody NBP2-47395 (A) shows similar protein distribution across tissues to independent antibody NBP2-47396 (B).



Immunohistochemistry-Paraffin: LMP7/PSMB8 Antibody [NBP2-47395] - Staining of human kidney.



Immunohistochemistry-Paraffin: LMP7/PSMB8 Antibody [NBP2-47395] - Staining of human liver.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-47395

NBP2-47395PEP	LMP7/PSMB8 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-47395

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

