

Product Datasheet

Cytokeratin, pan Antibody (PK110) [DyLight 488] NBP2-50222G

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-50222G

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-50222G



NBP2-50222G

Cytokeratin, pan Antibody (PK110) [DyLight 488]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	PK110
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 488
Purity	Protein A purified
Buffer	50mM Sodium Borate


Product Description	
Description	This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Host	Mouse
Gene ID	3848
Gene Symbol	KRT1
Species	Human
Specificity/Sensitivity	This antibody recognizes cytokeratins of 50-67 kDa.
Immunogen	This Cytokeratin, pan Antibody (PK110) was developed against epidermal keratin.
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Immunohistochemistry-Frozen
Application Notes	Optimal dilution of this antibody should be experimentally determined.




Images

Immunocytochemistry/Immunofluorescence: Cytokeratin, pan Antibody (PK110) [DyLight 488] [NBP2-50222G] - HaCaTs (keratinocytes) were stained with Cytokeratin, pan antibody (PK110) [DyLight 488]. Cells were fixed with PFA 4%. Antibody concentration: 1:100. Dilution buffer: 1X PBS, 1% bovine serum albumin (BSA). Image from verified customer review.



Immunohistochemistry-Frozen: Cytokeratin, pan Antibody (PK110) [DyLight 488] [NBP2-50222G] - Cytokeratin was detected in human skin using Cytokeratin, pan antibody (PK110) [DyLight 488] (1:100; 1 hour) in PBS. Tissue was fixed in acetone. Image from verified customer review.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-50222G

DNST0	Endostatin [HRP]
MAB1455	Albumin Antibody (188835) [Unconjugated] - Serum
NB100-687	Cytokeratin 19 Antibody - BSA Free
NBP2-16094	Cytokeratin 8 Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-50222G

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

