

Product Datasheet

RNA Polymerase II/POLR2A [p Ser5] Antibody (CTD 8A7) - BSA Free NBP2-50243

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-50243

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-50243



NBP2-50243

RNA Polymerase II/POLR2A [p Ser5] Antibody (CTD 8A7) - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CTD 8A7
Preservative	0.02% Sodium Azide
Isotype	IgM
Purity	Thiophilic Chromatography
Buffer	PBS


Product Description	
Description	Novus Biologicals Mouse RNA Polymerase II/POLR2A [p Ser5] Antibody (CTD 8A7) - BSA Free (NBP2-50243) is a monoclonal antibody validated for use in WB, ELISA, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5430
Gene Symbol	POLR2A
Species	Human, Yeast
Reactivity Notes	Saccharomyces cerevisiae
Immunogen	Peptide: 10 repeats of synthetic peptide YSPTSPS using chemically synthesized phospho-Ser5.

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, CyTOF-ready
Recommended Dilutions	Western Blot 1:100 - 1:2000, Flow Cytometry 1:10 - 1:1000, ELISA 1:100 - 1:2000, Immunocytochemistry/ Immunofluorescence 1:10 - 1:500, CyTOF-ready
Application Notes	This antibody is Cytof ready. Positive control(s): S. cerevisiae and HeLa cells.

Images

Western Blot: RNA Polymerase II/POLR2A Antibody (CTD 8A7) [NBP2-50243] - RNA polymerase II (RNAP II) contains 12 subunits. The C-terminal domain of the largest subunit of RNAP II is reversibly phosphorylated on the highly conserved heptapeptide repeat, YSPTSPS. The major physiological phosphorylation sites are Ser-2 and Ser-5. Monoclonal antibody 8A7 was obtained using synthetic phosphopeptides and specifically recognises the Ser-5 phosphorylated form of RNAP II.

Monoclonal antibody 8A7
 specific for hyperphosphorylated large subunit of RNA polymerase II
 cross-reactive: human, yeast
 mouse IgG
 use for immunoblotting, ELISA
 dilution for immunoblotting: 1:1,000





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-50243

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97007	Mouse IgM Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-50243

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

