

Product Datasheet

LONP1 Antibody (LON20/H1) - BSA Free NBP2-50478

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-50478

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-50478



NBP2-50478

LONP1 Antibody (LON20/H1) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	LON20/H1
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	PBS
Target Molecular Weight	106 kDa

Product Description	
Description	Novus Biologicals Mouse LONP1 Antibody (LON20/H1) - BSA Free (NBP2-50478) is a monoclonal antibody validated for use in IHC, WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	9361
Gene Symbol	LONP1
Species	Human
Immunogen	KLH-peptide conjugate - N-CEKDDKDAIEEKFRERLKE-C

Product Application Details	
Applications	Western Blot, ELISA, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot, ELISA, Immunohistochemistry, Knockdown Validated
Application Notes	Positive control(s): U87 cells

Images

Western Blot: LONP1 Antibody (LON20/H1) [NBP2-50478] - U87 cells were transfected with siRNA against the Lon protease (siLonp1.1 and siLonp1.2) or a control sequence (siControl) or transfection reagent alone (mock) and then western blotted using anti-LONP1 [LON20/H1] and a loading control antibody (B-actin). Molecular weight markers are labelled in kilodaltons





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-50478

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97019-5mg	Mouse IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-50478

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

