

# Product Datasheet

## CD30/TNFRSF8 Antibody (Ki-3) - Chimeric - Azide and BSA Free NBP2-52682-0.2mg

Unit Size: 0.2 mg

Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

### Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-52682](http://www.novusbio.com/NBP2-52682)

Updated 2/24/2026 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-52682](http://www.novusbio.com/reviews/destination/NBP2-52682)



**NBP2-52682-0.2mg**

CD30/TNFRSF8 Antibody (Ki-3) - Chimeric - Azide and BSA Free

Product Information	
Unit Size	0.2 mg
Concentration	1 mg/ml
Storage	Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.
Clonality	Monoclonal
Clone	Ki-3
Preservative	0.02% Proclin 300
Isotype	IgG Kappa
Purity	Protein A purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Rabbit CD30/TNFRSF8 Antibody (Ki-3) - Chimeric - Azide and BSA Free (NBP2-52682) is a recombinant monoclonal antibody validated for use in IHC, ELISA, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	943
Gene Symbol	TNFRSF8
Species	Human
Specificity/Sensitivity	The antibody binds to CD30.
Immunogen	Human CD30.

Product Application Details	
Applications	ELISA, Flow Cytometry, Functional, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Flow Cytometry 1:10 - 1:1000, ELISA 1:100 - 1:2000, Immunohistochemistry 1:10 - 1:500, Immunocytochemistry/ Immunofluorescence, Functional, CyTOF-ready
Application Notes	This antibody is Cytof ready. This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG2b format, for improved compatibility with existing reagents, assays and techniques. The antibody binds specifically to CD30/TNFRSF8 a cell membrane protein which belongs to the tumour necrosis factor receptor family. Incubation of CD30/TNFRSF8 positive L540 cells with Ki-3 antibody weakly enhances release of soluble CD30/TNFRSF8 (in contrast to Ki-4, NBP2-52714). CD30/TNFRSF8 is expressed on activated lymphocytes, and binding by its ligands TRAF2 and TRAF5 results in the activation of NF-kappaB, thereby regulating gene expression. CD30/TNFRSF8 is also associated with Hodgkin Lymphoma, and its expression acts as a tumour marker. Targeting the CD30/TNFRSF8 antigen with specific monoclonal antibodies is a possible therapy for Hodgkin's lymphoma.




## Images

**Immunohistochemistry: CD30/TNFRSF8 Antibody (Ki-3) - Chimeric [NBP2-52682]** - Immunohistochemical staining of human tonsil tissue using anti-CD30 antibody. Ki-3 Anti-CD30 staining of paraffin embedded human tonsil tissue using the rabbit-chimeric version of Ki-3 (NBP2-52682). Antigen retrieval was achieved by microwaving in citrate buffer, followed by blocking with protein block serum-free buffer. Primary antibody incubation with NBP2-52682 was carried out at 4 ug/ml for 30 minutes. Samples were then incubated with an anti-rabbit IgG HRP secondary antibody for 20 mins followed by DAB (3,3'-diaminobenzidine), and counter-staining with haematoxylin. Staining of the cellular membrane of dividing cells may be observed. Recommended concentration, 2-4 ug/ml.



**Flow Cytometry: CD30/TNFRSF8 Antibody (Ki-3) - Chimeric [NBP2-52682]** - Flow-cytometry using the anti-CD30 (Tumor necrosis factor receptor superfamily member 8) antibody Ki-3. Jurkat cells were stained with unimmunized rabbit IgG antibody (black line) or the rabbit-chimeric version of Ki-3 (NBP2-52682, blue line) at a concentration of 10 ug/ml for 30 mins at RT. After washing, bound antibody was detected using an anti-rabbit IgG JK (FITC-conjugate) antibody at 2 ug/ml and cells analyzed on a FACSCanto flow-cytometer.



## Publications

Sato K, Kozako T, Nakano A et al. Extracellular CD30 and ADAM10/17 regulate brentuximab vedotin-induced cell death in an adult T-cell leukemia cell line Research Square 2023-11-20 (ELISA, Human)



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-52682-0.2mg**

---

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-22660	Recombinant Human CD30/TNFRSF8 His Protein
210-TA-005	TNF-alpha [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-52682](http://www.novusbio.com/reviews/submit/NBP2-52682)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

