

Product Datasheet

TNF-alpha Antibody (TNF/1500R) NBP2-53242-100ug

Unit Size: 100 ug

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-53242

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-53242



NBP2-53242-100ug

TNF-alpha Antibody (TNF/1500R)

| Product Information | |
|-------------------------|--------------------------|
| Unit Size | 100 ug |
| Concentration | 0.2 mg/ml |
| Storage | Store at 4C. |
| Clonality | Monoclonal |
| Clone | TNF/1500R |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Purity | Protein A purified |
| Buffer | 10 mM PBS with 0.05% BSA |
| Target Molecular Weight | 17 kDa |


| Product Description | |
|---------------------|--|
| Description | 200ug/ml of antibody purified from Bioreactor Concentrate. Prepared in 10 mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0 mg/ml. (NBP3-11621) Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80C. |
| Host | Rabbit |
| Gene ID | 7124 |
| Gene Symbol | TNF |
| Species | Human, Rat |
| Immunogen | Recombinant full-length human TNF-alpha protein (Uniprot: P01375) |

| Product Application Details | |
|-----------------------------|---|
| Applications | Immunohistochemistry-Paraffin, Immunohistochemistry |
| Recommended Dilutions | Immunohistochemistry, Immunohistochemistry-Paraffin 4-8 ug/ml |
| Application Notes | Immunohistochemistry (Formalin-fixed): 4-8ug/ml for 30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes. Optimal dilution for a specific application should be determined. |




Images

Immunohistochemistry-Paraffin: TNF-alpha Antibody (TNF/1500R) [NBP2-53242] - Formalin-fixed, paraffin-embedded Human Pancreas stained with TNF-alpha Recombinant Rabbit Monoclonal Antibody (TNF/1500R).



Immunohistochemistry-Paraffin: TNF-alpha Antibody (TNF/1500R) [NBP2-53242] - Formalin-fixed, paraffin-embedded Human Histiocytoma stained with TNF-alpha Recombinant Rabbit Monoclonal Antibody (TNF/1500R).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-53242-100ug

| | |
|-----------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| NBP2-35076-10ug | Recombinant Human TNF-alpha Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-53242

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

