

# Product Datasheet

## CD43/Sialophorin Antibody (SPN/1766R) [mFluor Violet 500 SE] NBP2-54344MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-54344MFV500](http://www.novusbio.com/NBP2-54344MFV500)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-54344MFV500](http://www.novusbio.com/reviews/destination/NBP2-54344MFV500)



**NBP2-54344MFV500**

CD43/Sialophorin Antibody (SPN/1766R) [mFluor Violet 500 SE]

| <b>Product Information</b>         |   |
|------------------------------------|---|
| <b>Unit Size</b>                   | 0.1 ml  |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.   |
| <b>Storage</b>                     | Store at 4C in the dark.  |
| <b>Clonality</b>                   | Monoclonal  |
| <b>Clone</b>                       | SPN/1766R   |
| <b>Preservative</b>                | 0.05% Sodium Azide  |
| <b>Isotype</b>                     | IgG   |
| <b>Conjugate</b>                   | mFluor Violet 500 SE  |
| <b>Purity</b>                      | Protein A or G purified   |
| <b>Buffer</b>                      | 50mM Sodium Borate  |
| <b>Product Description</b>         |   |
| <b>Host</b>                        | Rabbit  |
| <b>Gene ID</b>                     | 6693  |
| <b>Gene Symbol</b>                 | SPN   |
| <b>Species</b>                     | Human   |
| <b>Marker</b>                      | T-Cell Marker   |
| <b>Specificity/Sensitivity</b>     | It recognizes a cell surface glycoprotein of 95/115/135kDa (depending upon the extent of glycosylation), identified as CD43. 70-90% of T-cell lymphomas and from 22-37% of B-cell lymphomas express CD43. No reactivity has been observed with reactive B-cells. So a B-lineage population that co-expresses CD43 is highly likely to be a malignant lymphoma, especially a low-grade lymphoma, rather than a reactive B-cell population. When CD43 antibody is used in combination with anti-CD20, effective immunophenotyping of the lymphomas in formalin-fixed tissues can be obtained. Co-staining of a lymphoid infiltrate with anti-CD20 and anti-CD43 argues against a reactive process and favors a diagnosis of lymphoma. |
| <b>Immunogen</b>                   | Recombinant full-length human CD43/Sialophorin protein (Uniprot: P16150)  |
| <b>Notes</b>                       | mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.   |
| <b>Product Application Details</b> |   |
| <b>Applications</b>                | Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, CyTOF-ready   |
| <b>Recommended Dilutions</b>       | Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, CyTOF-ready   |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.  |

## Images

CD43/Sialophorin Antibody (SPN/1766R) [mFluor Violet 500 SE] - Vial of mFluor Violet 500 conjugated antibody. mFluor Violet 500 is optimally excited at 410 nm by the Violet laser (405 nm) and has an emission maximum of 501 nm.



mFluor™ Violet 500

| LASER (nm)   | FILTER |
|--------------|--------|
| Violet (405) | 515/20 |

| EXCITATION MAX (nm) | EMISSION MAX (nm) |
|---------------------|-------------------|
| 410                 | 501               |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-54344MFV500**

---

|                  |   |
|------------------|---|
| NBP2-24891MFV500 | Rabbit IgG Isotype Control [mFluor Violet 500 SE] |
| NBP2-23472       | Recombinant Human CD43/Sialophorin His Protein    |
| 210-TA-005       | TNF-alpha [Unconjugated]                          |
| 9680-CD-050      | CD43/Sialophorin [Unconjugated]                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-54344MFV500](http://www.novusbio.com/reviews/submit/NBP2-54344MFV500)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

