

Product Datasheet

POU3F3 Antibody - BSA Free

NBP2-57011

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-57011

Updated 12/2/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-57011



NBP2-57011

POU3F3 Antibody - BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Rabbit POU3F3 Antibody - BSA Free (NBP2-57011) is a polyclonal antibody validated for use in IHC and ICC/IF. Anti-POU3F3 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5455
Gene Symbol	POU3F3
Species	Human, Mouse
Reactivity Notes	Use in Mouse reported in scientific literature (PMID:33649223).
Immunogen	This antibody was developed against a recombinant protein corresponding to the following amino acid sequence: GGAQSLVHPGLVRGDTPELAEHHHHH

Product Application Details	
Applications	Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Frozen
Recommended Dilutions	Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Frozen Reported in scientific literature (PMID:33649223)
Application Notes	ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.

Images

Immunocytochemistry/Immunofluorescence: POU3F3 Antibody [NBP2-57011] - Staining of human cell line SH-SY5Y shows localization to nucleus. Antibody staining is shown in green.



Publications

Tsujimoto, H;Hoshina, A;Mae, SI;Araoka, T;Changting, W;Ijiri, Y;Nakajima-Koyama, M;Sakurai, S;Okita, K;Mizuta, K;Niwa, A;Saito, MK;Saitou, M;Yamamoto, T;Graneli, C;Woollard, KJ;Osafune, K; Selective induction of human renal interstitial progenitor-like cell lineages from iPSCs reveals development of mesangial and EPO-producing cells Cell reports 2023-12-22 [PMID: 38237600]

Ruan X, Kang B, Qi C et al. Progenitor cell diversity in the developing mouse neocortex Proceedings of the National Academy of Sciences of the United States of America 2021-03-09 [PMID: 33649223] (IHC-Fr, Mouse)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-57011

NBP2-57011PEP	POU3F3 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-57011

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

