

Product Datasheet

TTC37 Antibody - BSA Free

NBP2-59165

Unit Size: 150 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-59165

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-59165



NBP2-59165

TTC37 Antibody - BSA Free

Product Information

Unit Size	150 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide and 0.05% ProClin 300
Isotype	IgG
Purity	Protein G purified
Buffer	PBS

Product Description

Description	Novus Biologicals Rabbit TTC37 Antibody - BSA Free (NBP2-59165) is a polyclonal antibody validated for use in WB, ICC/IF and ChIP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	9652
Gene Symbol	SKIC3
Species	Human
Immunogen	The exact immunogen is proprietary information.


Product Application Details

Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Chromatin Immunoprecipitation (ChIP), Immunofluorescence
Recommended Dilutions	Western Blot 1:1000, Immunocytochemistry/ Immunofluorescence 1:100, Immunofluorescence 1:100, Chromatin Immunoprecipitation (ChIP) 1-5 ug




Images

Western Blot: TTC37 Antibody [NBP2-59165] - Western blot was performed on HeLa nuclear extracts (lane 1) and on the purified human SKI complex (lane 2) (FigureA), or on nuclear extracts from HeLa cells (HeLa NE, 20 ug) (FigureB) using the antibody against SKI3 diluted 1:1,000. The location of the protein is indicated on the right, the marker (in kDa) is shown on the left.



Immunocytochemistry/Immunofluorescence: TTC37 Antibody [NBP2-59165] - HeLa cells were stained with the antibody against SKI3 (in red), with DAPI to label the DNA (in blue), and with an antibody against alpha-tubulin to label the cytoplasm (in green). SKI3 is detected in both the cytoplasm and the nucleus but appears to be excluded from the nucleolus.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-59165

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-59165

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

