

Product Datasheet

HA Tag Antibody (C.15200190) - BSA Free NBP2-59172

Unit Size: 50 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-59172

Updated 2/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-59172



NBP2-59172

HA Tag Antibody (C.15200190) - BSA Free


Product Information	
Unit Size	50 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	C.15200190
Preservative	0.05% Sodium Azide
Isotype	IgG3
Purity	Protein A purified
Buffer	PBS

Product Description	
Description	The HA antibody recognizes the influenza Hemagglutinin epitope which is extensively used as a protein tag in expression vectors. This antibody is extremely specific and allows unambiguous identification and quantitative analysis of HA tagged proteins.
Host	Mouse
Species	Epitope Tag
Immunogen	This HA Tag Antibody (C.15200190) was raised in mouse against the HA tag (amino acid sequence YPYDVPDYA) using a KLHconjugated synthetic peptide.

Product Application Details	
Applications	Western Blot, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP)
Recommended Dilutions	Western Blot 1:1000 - 1:3000, Immunoprecipitation 1:150 - 1:200, Chromatin Immunoprecipitation (ChIP) 1:150-1:200

Images

Western Blot: HA Tag Antibody (C.15200190) [NBP2-59172] - Whole cell extracts from HEK293 cells transfected with an HA-tagged expression vector using the antibody directed against HA at 1:1,000. Figure shows the results for untransfected cells, used as a negative control (lane 1) and for cells transfected with a HA-tagged protein (lane 2). The MW marker (in kDa) is shown on the left. A single band is clearly visible in lane 2, but absent in lane 1 showing the high specificity of the HA antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-59172

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96777	Mouse IgG3 Isotype Control (MG3)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-59172

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

