

# Product Datasheet

## LRRC35 Antibody (1B10) - BSA Free NBP2-59417-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-59417](http://www.novusbio.com/NBP2-59417)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-59417](http://www.novusbio.com/reviews/destination/NBP2-59417)




**NBP2-59417-100ul**

LRRC35 Antibody (1B10) - BSA Free


<b>Product Information</b>	
<b>Unit Size</b>	100 ul
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	1B10
<b>Preservative</b>	0.02% Sodium Azide
<b>Isotype</b>	IgG2b Kappa
<b>Purity</b>	Protein A purified
<b>Buffer</b>	PBS (pH 7.4), 10% Glycerol
<b>Product Description</b>	
<b>Description</b>	Novus Biologicals Mouse LRRC35 Antibody (1B10) - BSA Free (NBP2-59417) is a monoclonal antibody validated for use in WB, ELISA, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
<b>Host</b>	Mouse
<b>Gene ID</b>	219899
<b>Gene Symbol</b>	TBCEL
<b>Species</b>	Human
<b>Immunogen</b>	Recombinant human LRRC35 (1-424aa) purified from E. coli
<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence
<b>Recommended Dilutions</b>	Western Blot 1:1000, Flow Cytometry 1:10-1:1027, ELISA, Immunocytochemistry/ Immunofluorescence 1:100

## Images

Western Blot: LRRC35 Antibody (1B10) [NBP2-59417] - Lane 1: Recombinant Human TBCEL, Lane 2: PC3 cell lysates, Lane 3: LnCaP cell lysates, Lane 4: 293T cell lysates, Lane 5: Mouse testis lysates, Lane 6: Mouse Kidney lysates



Immunocytochemistry/Immunofluorescence: LRRC35 Antibody (1B10) [NBP2-59417] - Analysis of TBCEL in PC3 cells line, stained with DAPI (Blue) for nucleus staining and monoclonal anti-human TBCEL antibody (1:100) with goat anti-mouse IgG-Alexa fluor 488 conjugate (Green).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-59417-100ul**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-59417](http://www.novusbio.com/reviews/submit/NBP2-59417)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

