

Product Datasheet

Profilin 1 Antibody (CL3524) - BSA Free NBP2-59778

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-59778

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-59778



NBP2-59778

Profilin 1 Antibody (CL3524) - BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL3524
Preservative	0.02% Sodium Azide
Isotype	IgG2a
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Knockout (KO) Validated Mouse Profilin 1 Antibody (CL3524) - BSA Free (NBP2-59778) is a monoclonal antibody validated for use in IHC and WB. Anti-Profilin 1 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5216
Gene Symbol	PFN1
Species	Human
Reactivity Notes	Rat (87%), Mouse (87%).
Immunogen	This antibody was developed using a recombinant protein derived from P07737, with the exact immunogen sequence remaining proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Knockdown Validated, Knockout Validated
Recommended Dilutions	Western Blot 1 ug/ml, Immunohistochemistry 1:1000 - 1:2500, Immunohistochemistry-Paraffin 1:1000 - 1:2500, Knockout Validated, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Western Blot: Profilin 1 Antibody (CL3524) [NBP2-59778] - Profilin 1 Antibody (5216) [NBP2-59778] - Analysis of extracts from U-251 cells, transfected with: control siRNA, target specific siRNA probe #1, target specific siRNA probe #2, using Anti-PFN1 monoclonal antibody. Downregulation of antibody signal confirms target specificity. Remaining % intensity, relative control lane, is indicated. Anti-GAPDH monoclonal antibody was used as loading control.



Immunohistochemistry-Paraffin: Profilin 1 Antibody (5216) [NBP2-59778]
- Staining of human testis shows weak immunoreactivity in basal cells of the seminiferous ducts.




Western Blot: Profilin 1 Antibody (5216) [NBP2-59778] - Analysis in human cell line HeLa, human cell line HEK 293, human cell line A-431 and human cell line HepG2.



Immunohistochemistry-Paraffin: Profilin 1 Antibody (5216) [NBP2-59778]
- Staining of human fallopian tube shows cytoplasmic immunoreactivity in epithelial cells.



Immunohistochemistry-Paraffin: Profilin 1 Antibody (5216) [NBP2-59778]
- Staining of human tonsil shows moderate to strong immunoreactivity in lymphoid cells.




Immunohistochemistry-Paraffin: Profilin 1 Antibody (5216) [NBP2-59778] 


- Staining of human liver shows strong immunoreactivity in circulating lymphoid cells.

Immunohistochemistry-Paraffin: Profilin 1 Antibody (5216) [NBP2-59778]

- Staining of human kidney shows strong membranous combined with moderate cytoplasmic immunoreactivity in renal tubules.



Western Blot: Profilin 1 Antibody (CL3524) [NBP2-59778] - Lysates of HAP1 WT and PFN1 KO were prepared, and 50 ug of protein were processed for immunoblot with Profilin 1 antibody (NBP2-59778). The Ponceau stained transfer of the blot is shown. Antibody dilution used: NBP2-59778 at 1/1000. Predicted band size: 15 kDa. Image, protocol and testing courtesy of YCharOS Inc. (ycharos.com).



Publications

Ayoubi, R;McDowell, I;Fotouhi, M;Southern, K;McPherson, PS;Laflamme, C;NeuroSGC/YCharOS/EDDU collaborative group, ;ABIF Consortium, ; The identification of high-performing antibodies for Profilin-1 for use in Western blot, immunoprecipitation and immunofluorescence F1000Research 2023-03-29 [PMID: 37576538]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-59778

NBP2-59778PEP	Profilin 1 Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-59778

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

