

Product Datasheet

Human Plasminogen ELISA Kit (Colorimetric) NBP2-60515

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-60515

Updated 5/4/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-60515




NBP2-60515

Human Plasminogen ELISA Kit (Colorimetric)

| Product Information | |
|-----------------------------|--|
| Unit Size | 1 Kit |
| Concentration | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| Storage | Storage of components varies. See protocol for specific instructions. |
| Product Description | |
| Description | This assay employs a quantitative enzyme immunoassay technique that measures the specified antigen in samples. |
| Gene ID | 5340 |
| Gene Symbol | PLG |
| Species | Human |
| Standard Curve Range | 0.625 - 40 ng/mL (example only; lot dependent) |
| Sensitivity | 0.9 ng/ml (example only; lot dependent) |
| Inter Assay | CV% < 9.90% (example only; lot dependent) |
| Intra Assay | CV% < 5.00% (example only; lot dependent) |
| Assay Type | Sandwich ELISA |
| Spiking Recovery | 101.00% (example only; lot dependent) |
| Suitable Sample Type | Serum, Plasma, Cell Culture, Urine, Saliva, Milk, Tissue Extract, CSF |
| Product Application Details | |
| Applications | ELISA, Sandwich ELISA |
| Recommended Dilutions | ELISA, Sandwich ELISA |

Images

ELISA: Human Plasminogen ELISA Kit (Colorimetric) [NBP2-60515]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-60515

| | |
|-----------------|---|
| NBP1-86015PEP | Plasminogen Recombinant Protein Antigen |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| NB600-930-100ug | Plasminogen Antibody - BSA Free |
| M6000B-1 | IL-6 [HRP] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-60515

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

