

# Product Datasheet

## Hsp70 interacting protein HIP Antibody (CL3708) - BSA Free NBP2-61615

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-61615](http://www.novusbio.com/NBP2-61615)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-61615](http://www.novusbio.com/reviews/destination/NBP2-61615)



**NBP2-61615**

Hsp70 interacting protein HIP Antibody (CL3708) - BSA Free

<b>Product Information</b>	
<b>Unit Size</b>	100 ul
<b>Concentration</b>	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	CL3708
<b>Preservative</b>	0.02% Sodium Azide
<b>Isotype</b>	IgG1
<b>Purity</b>	Protein A purified
<b>Buffer</b>	PBS (pH 7.2) and 40% Glycerol


<b>Product Description</b>	
<b>Description</b>	Novus Biologicals Mouse Hsp70 interacting protein HIP Antibody (CL3708) - BSA Free (NBP2-61615) is a monoclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
<b>Host</b>	Mouse
<b>Gene ID</b>	6767
<b>Gene Symbol</b>	ST13
<b>Species</b>	Human
<b>Reactivity Notes</b>	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions
<b>Immunogen</b>	This antibody was developed using a recombinant protein derived from P50502, with the exact immunogen sequence remaining proprietary.

<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
<b>Recommended Dilutions</b>	Western Blot 1 ug/ml, Immunohistochemistry 1:10000 - 1:20000, Immunohistochemistry-Paraffin 1:10000 - 1:20000
<b>Application Notes</b>	IHC-Paraffin HIER pH6 retrieval is recommended.




## Images


Western Blot: Hsp70 interacting protein HIP Antibody (CL3708) [NBP2-61615] - Lane 1: Marker [kDa] 250, 130, 100, 70, 55, 35, 25, 15, 10  
Lane 2: Human cell line RT-4




Immunohistochemistry-Paraffin: Hsp70 interacting protein HIP Antibody (CL3708) [NBP2-61615] - Staining of human prostate cancer shows moderate cytoplasmic positivity in tumor cells.



Immunohistochemistry: Hsp70 interacting protein HIP Antibody (CL3708) [NBP2-61615] - Staining of human colorectal cancer shows cytoplasmic immunoreactivity in tumor cells.




Immunohistochemistry: Hsp70 interacting protein HIP Antibody (CL3708) [NBP2-61615] - Staining of human kidney shows cytoplasmic positivity in renal tubules and glomerulus cells.




Immunohistochemistry: Hsp70 interacting protein HIP Antibody (CL3708) [NBP2-61615] - Staining of human fallopian tube shows strong cytoplasmic immunoreactivity in epithelial cells.



Immunohistochemistry: Hsp70 interacting protein HIP Antibody (CL3708) [NBP2-61615] - Staining of human endometrium shows cytoplasmic immunoreactivity in glandular epithelium cells.



Immunohistochemistry-Paraffin: Hsp70 interacting protein HIP Antibody (CL3708) [NBP2-61615] - Staining of human liver shows cytoplasmic positivity in hepatocytes.



Immunohistochemistry-Paraffin: Hsp70 interacting protein HIP Antibody (CL3708) [NBP2-61615] - Staining of human breast cancer shows cytoplasmic immunoreactivity in tumor cells.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-61615**

---

NBP2-61615PEP	Hsp70 interacting protein HIP Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-61615](http://www.novusbio.com/reviews/submit/NBP2-61615)

Earn gift cards/discounts by submitting a publication using this product:

[www.novusbio.com/publications](http://www.novusbio.com/publications)

