

Product Datasheet

Fibronectin Antibody (CL3730) - BSA Free NBP2-61633

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-61633

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-61633



NBP2-61633**Fibronectin Antibody (CL3730) - BSA Free**

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL3730
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol


Product Description	
Description	Novus Biologicals Mouse Fibronectin Antibody (CL3730) - BSA Free (NBP2-61633) is a monoclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	2335
Gene Symbol	FN1
Species	Human
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions
Immunogen	This antibody was developed using a recombinant protein derived from P02751, with the exact immunogen sequence remaining proprietary.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1 ug/ml, Immunohistochemistry 1:1000 - 1:2500, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	IHC-Paraffin HIER pH6 retrieval is recommended.




Images


Western Blot: Fibronectin Antibody (CL3730) [NBP2-61633] - Lane 1: Marker [kDa] 250, 130, 100, 70, 55, 35, 25, 15, 10 Lane 2: Human cell line U-87 MG




Immunohistochemistry-Paraffin: Fibronectin Antibody (CL3730) [NBP2-61633] - Staining of human pancreas shows absence of immunoreactivity in glandular cells (negative control).




Immunohistochemistry-Paraffin: Fibronectin Antibody (CL3730) [NBP2-61633] - Staining of human endometrium shows strong immunoreactivity in endometrial stroma and basement membrane of glandular cells.




Immunohistochemistry-Paraffin: Fibronectin Antibody (CL3730) [NBP2-61633] - Staining of human placenta shows strong immunoreactivity in chorionic villi.




Immunohistochemistry-Paraffin: Fibronectin Antibody (CL3730) [NBP2-61633] - Staining of human testis shows immunoreactivity in the connective tissue surrounding seminiferous tubules.




Immunohistochemistry-Paraffin: Fibronectin Antibody (CL3730) [NBP2-61633] - Staining of human prostate shows moderate positivity in smooth muscle.




Immunohistochemistry-Paraffin: Fibronectin Antibody (CL3730) [NBP2-61633] - Staining of human small intestine shows strong immunoreactivity in the connective tissue.



Immunohistochemistry-Paraffin: Fibronectin Antibody (CL3730) [NBP2-61633] - Staining of human breast cancer shows strong immunoreactivity in tumor stroma.



Immunohistochemistry-Paraffin: Fibronectin Antibody (CL3730) [NBP2-61633] - Staining of human liver cancer shows strong positivity in tumor.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-61633

NBP2-61633PEP	Fibronectin Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-61633

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

