

Product Datasheet

AKAP12 Antibody - BSA Free NBP2-66604

Unit Size: 250 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-66604

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-66604



NBP2-66604

AKAP12 Antibody - BSA Free


| Product Information | |
|---------------------|---|
| Unit Size | 250 ug |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at -20C. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.08% Sodium Azide |
| Isotype | IgG |
| Purity | Ammonium sulfate precipitation |
| Buffer | PBS |

| Product Description | |
|---------------------|---|
| Description | Novus Biologicals Sheep AKAP12 Antibody - BSA Free (NBP2-66604) is a polyclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Sheep |
| Gene ID | 9590 |
| Gene Symbol | AKAP12 |
| Species | Rat |
| Immunogen | Sheep were immunized with a recombinant protein corresponding to amino acids 620-1596 of the SSeCKS protein (Uniprot: Q62766). |

| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot, ELISA, Immunoaffinity Purification |
| Recommended Dilutions | Western Blot 1:100-1:2000, ELISA, Immunoaffinity Purification |
| Application Notes | . |

Images

Western Blot: AKAP12 Antibody [NBP2-66604] - Western blotting of total cellular protein from cultured rat aortic smooth muscle cells was prepared and analyzed. The protein was transferred to a PVDF membrane, blocked with 1X TBS, 0.1% Tween 20/ 5% NFDM, probed with 5 ug/mL of pre-immune serum 5 ug of AKAP12 antibody (a 3 second exposure is shown). The blot was then washed and probed with monoclonal anti-sheep IgG coupled to HRP at 1:10000, and detected with luminol.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-66604

| | |
|-----------------|--|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP] |
| HAF016 | Donkey anti-Sheep IgG Secondary Antibody [HRP] |
| NL010 | Donkey anti-Sheep IgG Secondary Antibody [NL557] |
| NBP1-97055-10mg | Sheep IgG Isotype Control |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-66604

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

