

# Product Datasheet

## Haptoglobin Antibody (JM10-79) NBP2-66881

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-66881](http://www.novusbio.com/NBP2-66881)

Updated 2/24/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-66881](http://www.novusbio.com/reviews/destination/NBP2-66881)



**NBP2-66881****Haptoglobin Antibody (JM10-79)**


| Product Information  |  |
|----------------------|--|
| <b>Unit Size</b>     | 100 ul   |
| <b>Concentration</b> | 1 mg/ml  |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| <b>Clonality</b>     | Monoclonal   |
| <b>Clone</b>         | JM10-79  |
| <b>Preservative</b>  | 0.05% Sodium Azide   |
| <b>Isotype</b>       | IgG  |
| <b>Purity</b>        | Protein A purified   |
| <b>Buffer</b>        | TBS (pH7.4), 0.05% BSA, 40% Glycerol   |

| Product Description |  |
|---------------------|--|
| <b>Description</b>  | Novus Biologicals Rabbit Haptoglobin Antibody (JM10-79) (NBP2-66881) is a recombinant monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| <b>Host</b>         | Rabbit   |
| <b>Gene ID</b>      | 3240   |
| <b>Gene Symbol</b>  | HP   |
| <b>Species</b>      | Human, Mouse, Rat  |
| <b>Immunogen</b>    | Synthetic peptide within Human Haptoglobin aa 357-406 / 406. (SwissProt: P00738 Human)   |


| Product Application Details  |  |
|------------------------------|--|
| <b>Applications</b>          | Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry                                    |
| <b>Recommended Dilutions</b> | Western Blot 1:500-1:2000, Immunohistochemistry, Immunocytochemistry/Immunofluorescence 1:50-1:200, Immunohistochemistry-Paraffin 1:50-1:200 |

**Images**

Western Blot: Haptoglobin Antibody (JM10-79) [NBP2-66881] - Analysis of Haptoglobin on different cells lysates using anti-Haptoglobin antibody at 1/500 dilution. Positive control: Lane 1: Hela Lane 2: HepG2




Immunocytochemistry/Immunofluorescence: Haptoglobin Antibody (JM10-79) [NBP2-66881] - Staining Haptoglobin in HepG2 cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunohistochemistry-Paraffin: Haptoglobin Antibody (JM10-79) [NBP2-66881] - Analysis of paraffin-embedded mouse lung tissue using anti-Haptoglobin antibody. Counter stained with hematoxylin.




Immunocytochemistry/Immunofluorescence: Haptoglobin Antibody (JM10-79) [NBP2-66881] - Staining Haptoglobin in A549 cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunocytochemistry/Immunofluorescence: Haptoglobin Antibody (JM10-79) [NBP2-66881] - Staining Haptoglobin in Hela cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunohistochemistry-Paraffin: Haptoglobin Antibody (JM10-79) [NBP2-66881] - Analysis of paraffin-embedded human liver tissue using anti-Haptoglobin antibody. Counter stained with hematoxylin.



Immunohistochemistry-Paraffin: Haptoglobin Antibody (JM10-79) [NBP2-66881] - Analysis of paraffin-embedded human lung tissue using anti-Haptoglobin antibody. Counter stained with hematoxylin.



Immunohistochemistry-Paraffin: Haptoglobin Antibody (JM10-79) [NBP2-66881] - Analysis of paraffin-embedded mouse liver tissue using anti-Haptoglobin antibody. Counter stained with hematoxylin.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-66881**

---

|             |   |
|-------------|---|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP]                 |
| HAF008      | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160      | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891  | Rabbit IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-66881](http://www.novusbio.com/reviews/submit/NBP2-66881)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

