

Product Datasheet

MUC2 Antibody (JA10-92)

NBP2-66961

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-66961

Updated 2/24/2026 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-66961



NBP2-66961

MUC2 Antibody (JA10-92)


Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	JA10-92
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	TBS (pH7.4), 0.05% BSA, 40% Glycerol
Target Molecular Weight	540 kDa

Product Description	
Description	Novus Biologicals Rabbit MUC2 Antibody (JA10-92) (NBP2-66961) is a recombinant monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. Anti-MUC2 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	4583
Gene Symbol	MUC2
Species	Human
Specificity/Sensitivity	MUC2 Antibody (JA10-91) is tissue specific for colon, small intestine, colonic tumors, bronchus, cervix and gall bladder.
Immunogen	Synthetic peptide within Human MUC2 aa 5000-5049 / 5179. (SwissProt: Q02817 Human; SwissProt: Q62635 Rat)


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:1000, Flow Cytometry 1:50-1:100, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:100-1:500, Immunohistochemistry-Paraffin 1:50-1:200

Images


Western Blot: MUC2 Antibody (JA10-92) [NBP2-66961] - Analysis of MUC2 on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1:5,000 dilution was used for 1 hour at room temperature. Positive control: Lane 1: SHSY-5Y cell lysate Lane 2: SKBR3 cell lysate




Immunocytochemistry/Immunofluorescence: MUC2 Antibody (JA10-92) [NBP2-66961] - Staining MUC2 in SW480 cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunohistochemistry-Paraffin: MUC2 Antibody (JA10-92) [NBP2-66961] - Immunohistochemical analysis of paraffin-embedded human small intestine tissue using anti-MUC2 antibody. The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 8.0-8.4) for 20 minutes. The tissues were blocked in 5% BSA for 30 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (1/200) for 30 minutes at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.




Flow Cytometry: MUC2 Antibody (JA10-92) [NBP2-66961] - Analysis of HeLa cells with MUC2 antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti rabbit IgG was used as the secondary antibody.




Immunocytochemistry/Immunofluorescence: MUC2 Antibody (JA10-92) [NBP2-66961] - Staining MUC2 in HeLa cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunocytochemistry/Immunofluorescence: MUC2 Antibody (JA10-92) [NBP2-66961] - Staining MUC2 in HepG2 cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunohistochemistry-Paraffin: MUC2 Antibody (JA10-92) [NBP2-66961] - Analysis of paraffin-embedded human breast cancer tissue using anti-MUC2 antibody. Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: MUC2 Antibody (JA10-92) [NBP2-66961] - Analysis of paraffin-embedded human colon cancer tissue using anti-MUC2 antibody. Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: MUC2 Antibody (JA10-92) [NBP2-66961] - Analysis of paraffin-embedded human tonsil tissue using anti-MUC2 antibody. Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: MUC2 Antibody (JA10-92) [NBP2-66961] - Analysis of paraffin-embedded rat brain tissue using anti-MUC2 antibody. Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: MUC2 Antibody (JA10-92) [NBP2-66961] - Analysis of paraffin-embedded rat intestine tissue using anti-MUC2 antibody. Counter stained with hematoxylin.




Immunohistochemistry-Paraffin: MUC2 Antibody (JA10-92) [NBP2-66961] - Analysis of paraffin-embedded rat skeletal muscle tissue using anti-MUC2 antibody. Counter stained with hematoxylin.



Immunohistochemistry-Paraffin: MUC2 Antibody (JA10-92) [NBP2-66961] - Immunofluorescence staining of paraffin- embedded human colon using anti-MUC2 rabbit polyclonal antibody. The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 10% negative goat serum for 1 hour at room temperature, washed with PBS, and then probed with at 1/100 dilution for 10 hours at 4C and detected using Alexa Fluor(R) 488 conjugate-Goat anti-Rabbit IgG (H+L) Secondary Antibody at a dilution of 1:500 for 1 hour at room temperature.



Immunohistochemistry-Paraffin: MUC2 Antibody (JA10-92) [NBP2-66961] - Immunohistochemical analysis of paraffin-embedded human colon tissue using anti-MUC2 antibody. The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 8.0-8.4) for 20 minutes. The tissues were blocked in 5% BSA for 30 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (1/200) for 30 minutes at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.



Publications

Fois CAM, Schindeler A, Valtchev P, Dehghani F Dynamic flow and shear stress as key parameters for intestinal cells morphology and polarization in an organ-on-a-chip model Biomedical microdevices 2021-10-16 [PMID: 34655329] (ICC/IF)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-66961

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-66961

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

