

Product Datasheet

alpha-Fetoprotein/AFP Antibody (OTI2A9) - Azide and BSA Free NBP2-70103

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-70103

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-70103




NBP2-70103

alpha-Fetoprotein/AFP Antibody (OTI2A9) - Azide and BSA Free


Product Information	
Unit Size	100 ug
Concentration	LYOPH mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI2A9
Preservative	No Preservative
Reconstitution Instructions	we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Purity	Immunogen affinity purified
Buffer	Lyophilized from PBS (pH 7.3) with 8% Trehalose
Target Molecular Weight	66.3 kDa
Product Description	
Description	Novus Biologicals Mouse alpha-Fetoprotein/AFP Antibody (OTI2A9) - Azide and BSA Free (NBP2-03084) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	174
Gene Symbol	AFP
Species	Human, Canine, Primate, Monkey
Immunogen	Full length human recombinant protein of human AFP (NP_001125) produced in HEK293T cell.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Western Blot 1:200-500, Flow Cytometry 1:100, Immunohistochemistry 1:150, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin, CyTOF-ready

Images

Western Blot: alpha-Fetoprotein/AFP Antibody (OT12A9) - Azide and BSA Free [NBP2-70103] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY Alpha 1 Fetoprotein (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-A1




Immunocytochemistry/Immunofluorescence: alpha-Fetoprotein/AFP Antibody (OT12A9) - Azide and BSA Free [NBP2-70103] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY Alpha 1 Fetoprotein.




Immunohistochemistry: alpha-Fetoprotein/AFP Antibody (OT12A9) - Azide and BSA Free [NBP2-70103] - Staining of paraffin-embedded Human liver tissue using anti-Alpha 1 Fetoprotein mouse monoclonal antibody.




Flow Cytometry: alpha-Fetoprotein/AFP Antibody (OT12A9) - Azide and BSA Free [NBP2-70103] - HEK293T cells transfected with either overexpression plasmid (Red) or empty vector control plasmid (Blue) were immunostaining by anti-Alpha 1 Fetoprotein antibody, and then analyzed by flow cytometry.



Western Blot: alpha-Fetoprotein/AFP Antibody (OT12A9) - Azide and BSA Free [NBP2-70103] - Analysis of extracts (35ug) from 9 different cell lines by using anti-Alpha 1 Fetoprotein monoclonal antibody.



Flow Cytometry: alpha-Fetoprotein/AFP Antibody (OT12A9) - Azide and BSA Free [NBP2-70103] - Analysis of Jurkat cells, using anti-Alpha 1 Fetoprotein antibody, (Red), compared to a nonspecific negative control antibody (Blue).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-70103

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-70103

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

