

Product Datasheet

Eukaryotic translation initiation factor 1 Antibody (OTI2B9) - Azide and BSA Free NBP2-71554

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-71554

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-71554




NBP2-71554**Eukaryotic translation initiation factor 1 Antibody (OTI2B9) - Azide and BSA Free**


Product Information	
Unit Size	100 ug
Concentration	LYOPH mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI2B9
Preservative	No Preservative
Reconstitution Instructions	we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Purity	Immunogen affinity purified
Buffer	Lyophilized from PBS (pH 7.3) with 8% Trehalose
Target Molecular Weight	12.6 kDa
Product Description	
Description	Novus Biologicals Mouse Eukaryotic translation initiation factor 1 Antibody (OTI2B9) - Azide and BSA Free (NBP2-00702) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10209
Gene Symbol	EIF1
Species	Human, Rat
Immunogen	Full length human recombinant protein of human EIF1 (NP_005792) produced in HEK293T cell.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Western Blot 1:2000, Flow Cytometry 1:100, Immunohistochemistry 1:150, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin, CyTOF-ready

Images


Western Blot: Eukaryotic translation initiation factor 1 Antibody (OTI2B9) - Azide and BSA Free [NBP2-71554] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY Eukaryotic translation initiation factor 1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and imm




Immunocytochemistry/Immunofluorescence: Eukaryotic translation initiation factor 1 Antibody (OTI2B9) - Azide and BSA Free [NBP2-71554] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY Eukaryotic translation initiation factor 1.




Immunohistochemistry: Eukaryotic translation initiation factor 1 Antibody (OTI2B9) - Azide and BSA Free [NBP2-71554] - Staining of paraffin-embedded Human lymph node tissue using anti-Eukaryotic translation initiation factor 1 mouse monoclonal antibody.




Flow Cytometry: Eukaryotic translation initiation factor 1 Antibody (OTI2B9) - Azide and BSA Free [NBP2-71554] - Analysis of Jurkat cells, using anti-Eukaryotic translation initiation factor 1 antibody, (Red), compared to a nonspecific negative control antibody (Blue).




Immunocytochemistry/Immunofluorescence: Eukaryotic translation initiation factor 1 Antibody (OTI2B9) - Azide and BSA Free [NBP2-71554] - Staining of HeLa cells using anti-Eukaryotic translation initiation factor 1 mouse monoclonal antibody.



Flow Cytometry: Eukaryotic translation initiation factor 1 Antibody (OTI2B9) - Azide and BSA Free [NBP2-71554] - Analysis of HeLa cells, using anti-Eukaryotic translation initiation factor 1 antibody, (Red), compared to a nonspecific negative control antibody (Blue).



CytoTOF-ready: Eukaryotic translation initiation factor 1 Antibody (OTI2B9) - Azide and BSA Free [NBP2-71554] - HEK293T cells transfected with either overexpression plasmid (Red) or empty vector control plasmid (Blue) were immunostaining by anti-Eukaryotic translation initiation factor 1 antibody, and then analyzed by flow cytometry.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-71554

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-71554

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

