

Product Datasheet

Indoleamine 2,3-dioxygenase/IDO Antibody (BLR040F) NBP2-76408

Unit Size: 100 ul

Store at 2 - 8 C / 1 year from date of receipt

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-76408

Updated 2/24/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-76408




NBP2-76408**Indoleamine 2,3-dioxygenase/IDO Antibody (BLR040F)**


Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 2 - 8 C / 1 year from date of receipt
Clonality	Monoclonal
Clone	BLR040F
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	Borate Buffered Saline (BBS) pH 8.2 with 0.1% BSA
Product Description	
Description	Novus Biologicals Rabbit Indoleamine 2,3-dioxygenase/IDO Antibody (BLR040F) (NBP2-76408) is a recombinant monoclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3620
Gene Symbol	IDO1
Species	Human
Immunogen	Recombinant protein representing Human IDO1
Notes	This Product sold under License from Bethyl Laboratories, Inc.
Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000, Immunoprecipitation 15 ul/1mg lysate

Images

Western Blot: Indoleamine 2,3-dioxygenase/IDO Antibody (BLR040F) [NBP2-76408] - Detection of human IDO1 by WB of HeLa treated with IFN-gamma (+) or mock treated (-). Antibody: Rabbit anti-IDO1 recombinant monoclonal [BLR040F] (NBP2-76408). Secondary: HRP-conjugated goat anti-rabbit IgG. Lower Panel: Rabbit anti-Vinculin.



Western Blot: Indoleamine 2,3-dioxygenase/IDO Antibody (BLR040F) [NBP2-76408] - Detection of human IDO1 by WB of immunoprecipitates from HeLa lysate treated with IFN-gamma. Antibodies: Rabbit anti-IDO1 recombinant monoclonal [BLR040F] (NBP2-76408) and rabbit anti-IDO1 recombinant monoclonal [BLR031F] (NBP2-76404).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-76408

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-76408

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

