

Product Datasheet

CD155/PVR Antibody (BLR074G) NBP2-78711

Unit Size: 100 uL

Store at 2 - 8 C / 1 year from date of receipt

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-78711

Updated 2/24/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-78711



NBP2-78711


CD155/PVR Antibody (BLR074G)

| Product Information | |
|------------------------------------|--|
| Unit Size | 100 uL |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 2 - 8 C / 1 year from date of receipt |
| Clonality | Monoclonal |
| Clone | BLR074G |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Protein A or G purified from cell culture supernatant |
| Buffer | Borate Buffered Saline (BBS) pH 8.2 with 0.1% BSA. |
| Product Description | |
| Description | Novus Biologicals Rabbit CD155/PVR Antibody (BLR074G) (NBP2-78711) is a recombinant monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Rabbit |
| Gene ID | 5817 |
| Gene Symbol | PVR |
| Species | Human |
| Immunogen | A recombinant protein representing the extracellular domain of human PVR/CD155. |
| Notes | This Product sold under License from Bethyl Laboratories, Inc. |
| Product Application Details | |
| Applications | Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry |
| Recommended Dilutions | Western Blot 1:1000, Flow Cytometry Fixed in 4% formaldehyde and permeabilized with 90% methanol. 2 µl, Immunohistochemistry 1:100 - 1:500, Immunocytochemistry/ Immunofluorescence 1:100 - 1:500, Immunohistochemistry-Paraffin 1:100 - 1:500 |
| Application Notes | Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use. IHC and ICC/IF: Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE cell sections. |




Images


Western Blot: CD155/PVR Antibody (BLR074G) [NBP2-78711] - Detection of human PVR/CD155 by WB of 786-O, K-562, OVCAR-4, Jurkat, KM12, NCI-H226, RPMI-8226, and NCC-IT lysate using this antibody. Secondary: HRP-conjugated goat anti-rabbit IgG (NB7160). Lower Panel: Rabbit anti-GAPDH (NB300-322).




Immunocytochemistry/Immunofluorescence: CD155/PVR Antibody (BLR074G) [NBP2-78711] - Detection of human PVR/CD155 in FFPE HT29 cells by ICC using this antibody. Secondary: HRP-conjugated goat anti-rabbit IgG (NB100-2332). Substrate: DAB.




Immunohistochemistry-Paraffin: CD155/PVR Antibody (BLR074G) [NBP2-78711] - Antibody: Rabbit anti-PVR/CD155 recombinant monoclonal antibody [BLR074G]. Secondary: HRP-conjugated goat anti-rabbit IgG Substrate: Opal™. Counterstain: DAPI.



Immunohistochemistry-Paraffin: CD155/PVR Antibody (BLR074G) [NBP2-78711] - Detection of human PVR/CD155 in FFPE lung carcinoma by IHC using this antibody. Secondary: HRP-conjugated goat anti-rabbit IgG (NB100-2332). Substrate: DAB.



Flow Cytometry: CD155/PVR Antibody (BLR074G) [NBP2-78711] -: Rabbit anti-PVR/CD155 recombinant monoclonal antibody [BLR074G] or isotype control (unshaded). Secondary: DyLight® 488- conjugated goat anti-rabbit IgG.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-78711

| | |
|-------------|---|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP] |
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-78711

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

