

Product Datasheet

GITR Ligand/TNFSF18 Antibody - Azide and BSA Free NBP2-80104

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80104

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80104



NBP2-80104**GITR Ligand/TNFSF18 Antibody - Azide and BSA Free**

Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	0.2 um-filtered solution in PBS, pH 7.4


Product Description	
Description	Novus Biologicals Rabbit GITR Ligand/TNFSF18 Antibody - Azide and BSA Free (NBP2-80104) is a polyclonal antibody validated for use in IHC, WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	8995
Gene Symbol	TNFSF18
Species	Human
Specificity/Sensitivity	Recognizes human GITR Ligand/TNFSF18. Detects a band of ~15-17kDa by Western blot. Does not cross-react with mouse GITR Ligand/TNFSF18.
Immunogen	Recombinant human GITR Ligand/TNFSF18.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry
Recommended Dilutions	Western Blot 1:2000 - 1:5000, ELISA 1:5000 - 1:10000, Immunohistochemistry 5 ug/mL, Immunohistochemistry-Paraffin 5 ug/mL




Images

Immunohistochemistry-Paraffin: GITR Ligand/TNFSF18 Antibody [NBP2-80104] - Immunohistochemical staining of GITRL using anti-GITRL (human), pAb in human uterus tissue (5 ug/mL). This antibody has been tested in immunohistochemistry, analyzed by an anatomic pathologist and validated for use in IHC applications against formalin-fixed, paraffin-embedded human tissues. The image shows the localization of the antibody as the precipitated red signal, with a hematoxylin purple nuclear counterstain.



GITR Ligand/TNFSF18 Antibody [NBP2-80104] - Western blot analysis using anti-GITRL (human), pAb at 1:5000 dilution. 1: Human GITRL (His-tagged). 2: Human GITRL (FLAG(R)-tagged). 3: Mouse GITRL (FLAG(R)-tagged).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80104

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80104

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

