

Product Datasheet

3-OHKYN Antibody (P3UI) - Azide and BSA Free NBP2-80548

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80548

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80548



NBP2-80548

3-OHKYN Antibody (P3UI) - Azide and BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	2.6 mg/ml
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	P3UI
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Protein G purified
Buffer	Tris-Glycine, 0.15 M NaCl
Target Molecular Weight	66 kDa


Product Description	
Description	Novus Biologicals Mouse 3-OHKYN Antibody (P3UI) - Azide and BSA Free (NB100-597) is a monoclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Species	Human, All Species
Immunogen	3-OHKYN modified KLH

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 2 ug/ml, ELISA reported in scientific literature, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:100-1:200, Immunohistochemistry-Paraffin reported in scientific literature




Images

Western Blot: 3-OHKYN Antibody (P3UI) - Azide and BSA Free [NBP2-80548] - Detection of 3'OHKYN in modified BSA using NB 100-597. 10 second ECL exposure. Lane 1: BSA (-) Lane 2: BSA (+). Image from the standard format of this antibody.



Immunohistochemistry-Paraffin: 3-OHKYN Antibody (P3UI) - Azide and BSA Free [NBP2-80548] - 3-OHKYN was detected in immersion fixed paraffin-embedded sections of human brain using Mouse Anti-Human 3-OHKYN (P3UI) Monoclonal Antibody (Catalog # NB100-597) at 1:300 for 1 hour at room temperature followed by incubation with the Anti-Mouse IgG VisUCy





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80548

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80548

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

